

## Abrasive Mounted Stones

### FG mounted stones



### FG GREEN ARKANSAS MOUNTED STONES - BESQUAL

#### FG GREEN ARKANSAS MOUNTED STONES - BESQUAL

Silicon Carbide For Porcelain, Gold, Amalgam and Enamel



Code	Name	Price	+
101770	BESQUAL FG CN1 GREEN Mounted Stones - Cone Shape 12pk - Silicone carbide for por...	US\$ 8.69	<a href="#">+Info</a>
102224	BESQUAL FG CY1 GREEN Mounted Stones - Small Cylinder Shape 12pk - Silicone carbi...	US\$ 7.90	<a href="#">+Info</a>
102225	BESQUAL FG CY2 GREEN Mounted Stones - Large Cylinder Shape 12pk - Silicone carbi...	US\$ 8.69	<a href="#">+Info</a>
102226	BESQUAL FG FL2 GREEN Mounted Stones - Flame Shape 12pk - Silicone carbide for po...	US\$ 9.87	<a href="#">+Info</a>
102227	BESQUAL FG IC1 GREEN Mounted Stones - Inverted Cone Shape 12pk - Silicone carbi...	US\$ 8.39	<a href="#">+Info</a>
102228	BESQUAL FG TC1 GREEN Mounted Stones - Taper Shape 12pk - Silicone carbide for ...	US\$ 9.87	<a href="#">+Info</a>

### FG WHITE MOUNTED STONES - BESQUAL

#### FG WHITE MOUNTED STONES - BESQUAL

For fine finishing of composites



#### BESQUAL

Code	Name	Price	+
102233	FG White Mnt Stone In Con 12pk	US\$ 7.90	<a href="#">+Info</a>

#### JSP












Code	Name	Price	+
101771	FG White Mnt Stone Cone 12pk	US\$ 7.15	<a href="#">+Info</a>
102229	FG White Mnt Stone Cyl 12pk	US\$ 7.01	<a href="#">+Info</a>
102230	FG White Mnt Stone Flame 12pk	US\$ 7.15	<a href="#">+Info</a>
102231	FG White Mnt Stone LFlame 12pk	US\$ 7.15	<a href="#">+Info</a>
102232	FG White Mnt Stone Rnd 12pk	US\$ 7.15	<a href="#">+Info</a>

## ABRASIVE MOUNTED STONES - CARGUS



Code	Name	Price	+
425520	MARK3 Cargus FL2 flame FG (friction)	US\$ 11.81	<a href="#">+Info</a>
425553	MARK3 Cargus CN1 pointed cone FG (f	US\$ 11.81	<a href="#">+Info</a>
425604	MARK3 Cargus RD1 round FG (friction)	US\$ 11.81	<a href="#">+Info</a>

## FG Friction Grip

SHAPE CATEGORY	CN1	CN2	CN3	CY2	CY3	CY4	FL1	FL2	FL3	FL4	FL5
SHAPE											
ISO #	024	037	040	025	028	045	030	025	030	045	065
CA (12 pk)	PN 0052	PN 0053	PN 0054	-	PN 0057	-	PN 0059	PN 0060	PN 0061	PN 0062	-
FG (12 pk)	PN 0101	-	-	PN 0103	-	-	-	PN 0104	-	-	-

## Cylindrical

Code	Name	Price	+
111485	SHOFU FG DURA-GREEN Cylindrical Cy2 12pk - #0103	US\$ 21.64	<a href="#">+Info</a>












## Flame

Code	Name	Price	+
111492	SHOFU Fg DURA-GREEN Flame FI1 12pk - #0009	US\$ 79.66	<a href="#">+Info</a>
111486	SHOFU FG DURA-GREEN Flame FI2 12pk - #0104	US\$ 21.64	<a href="#">+Info</a>

## Inverted Cone

Code	Name	Price	+
111487	SHOFU FG DURA-GREEN Inverted Cone IC2 12pk - #0105	US\$ 21.64	<a href="#">+Info</a>

## Round

SHAPE CATEGORY	KN4	KN5	KN6	KN7	PC1	PC2	RD1	RD2	RD3	RD4	RE1
SHAPE											
ISO #	080	060	100	125	030	050	030	024	036	050	062
CA (12 pk)	PN -	PN -	PN -	PN -	PN 0080	PN 0081	PN 0082	PN 0083	PN 0084	PN 0085	PN 0086
FG (12 pk)	PN -	PN -	PN -	PN -	PN -	PN -	PN 0107	PN 0108	PN -	PN -	PN -

Code	Name	Price	+
111488	SHOFU FG DURA-GREEN Round Rd1 12pk - #0107	US\$ 21.64	<a href="#">+Info</a>
111489	SHOFU FG DURA-GREEN Round Rd2 12pk - #0108	US\$ 21.21	<a href="#">+Info</a>


## Taper Conic

Code	Name	Price	+
111490	SHOFU FG DURA-GREEN Tc1 12pk - #0109	US\$ 21.64	<a href="#">+Info</a>

## Wheels

Code	Name	Price	+
111491	SHOFU FG DURA-GREEN Wheel WH1 12pk - #0110	US\$ 21.64	<a href="#">+Info</a>

## FG Friction Grip Stones

SHAPE CATEGORY	WH1
SHAPE	
ISO #	040
FG (12 PK)	PN 0250

Code	Name	Price	+
111580	SHOFU FG Dura-WHITE Wheel WH1 12pk - #0250	US\$ 22.30	<a href="#">+Info</a>
111579	SHOFU FG Dura-WHITE Taper Conic 1 12pk - #0249	US\$ 22.30	<a href="#">+Info</a>
111578	SHOFU FG Dura-WHITE Round Rd2 12pk - #0248	US\$ 19.85	<a href="#">+Info</a>
111577	SHOFU FG Dura-WHITE Round Rd1 12pk - #0247	US\$ 21.85	<a href="#">+Info</a>
111576	SHOFU FG Dura-WHITE Knife Kn1 12pk - #0246	US\$ 22.30	<a href="#">+Info</a>
111575	SHOFU FG Dura-WHITE Inverted Cone IC2 12pk - #0245	US\$ 22.30	<a href="#">+Info</a>
111574	SHOFU FG Dura-WHITE Flame FI2 12pk - #0244	US\$ 22.08	<a href="#">+Info</a>
111573	SHOFU FG Dura-WHITE Cylindrical Cy2 12pk - #0243	US\$ 22.08	<a href="#">+Info</a>
111572	SHOFU FG Dura-WHITE Cone Cn112pk - #0241	US\$ 19.85	<a href="#">+Info</a>
111571	SHOFU FG Dura-WHITE Flame FL1 12pk - #0240	US\$ 20.52	<a href="#">+Info</a>

## FG mounted stones - NEW BRAND SAMPLE

FG mounted stones - NEW BRAND SAMPLE

### Abrasive Mounted Stones

## Green mounted points



## DURA GREEN CA - SHOFU

### DURA GREEN - FOR CLINIC CA - SHOFU

#### Advantage:

Dura-Green Stones are excellent for fast contouring and finishing. Dura-Green Fine Stones are ideal for surface characterization and adjusting details.

Impregnated with silicon carbide. Long lasting, fast cutting and always concentric in motion.

Use Dura-Green Stones for: composites, amalgams, precious alloys, non-precious alloys and porcelains.

Use Dura-Green Fine Stones for: composites and porcelains See page 18, Unmounted Finishing Wheels section, for Dura-Green Wheels.














### Shofu raves and recommendations:

Sharpening Dental Instruments with Dura-Green? or Dura-White? Stones

After years of hit and miss sharpening of dental hygiene instruments, Dianne Glasscoe Watterson has developed a technique that is easier, faster, more effective, and promotes better ergonomics than any other manual technique. Using a high-speed handpiece, Shofu Dura-Green or a Dura-White stones, plus magnification, Dianne Glasscoe Watterson visually examines the instrument and removes the shiny, dull edge part of the blade until the cutting edge does not reflect light.

US DENTAL DEPOT Inc is an Authorized SHOFU Dealer in the US.

### Cone

SHAPE CATEGORY	CN1	CN2	CN3	CY2	CY3	CY4	FL1	FL2	FL3	FL4	FL5
SHAPE											
ISO #	024	037	040	025	028	045	030	025	030	045	065
CA (12 pk)	PN 0052	PN 0053	PN 0054	-	PN 0057	-	PN 0059	PN 0060	PN 0061	PN 0062	-
FG (12 pk)	PN 0101	-	-	PN 0103	-	-	-	PN 0104	-	-	-

Code	Name	Price	+
111462	SHOFU CA DURA-GREEN Cone Cn1 12pk - #0052	US\$ 19.48	<a href="#">+Info</a>
111463	SHOFU CA DURA-GREEN Cone Cn2 12pk - #0053	US\$ 21.64	<a href="#">+Info</a>
111464	SHOFU CA DURA-GREEN Cone Cn1 12pk - #0054	US\$ 21.64	<a href="#">+Info</a>












### Cylindrical

Code	Name	Price	+
111465	SHOFU CA DURA-GREEN Cylindrical Cy3 12pk - #0057	US\$ 21.64	<a href="#">+Info</a>

### Flame












Code	Name	Price	+
111466	SHOFU CA DURA-GREEN Flame Fl3 12pk - #0061	US\$ 21.64	<a href="#">+Info</a>

### Inverted Cone

SHAPE CATEGORY	IC1	IC2	IC3	IC4	IC5	IC6	IC7	IC8	IC9	KN1	KN3
SHAPE											
ISO #	030	032	050	045	070	065	125	104	135	050	070
CA (12 pk)	PN 0064	-	PN 0066	-	PN 0068	-	-	-	-	-	-
FG (12 pk)	PN -	PN 0105	-	-	-	-	-	-	-	-	-

Code	Name	Price	+
111428	SHOFU Ca DURA-GREEN Inverted Cone IC2 12pk - #0015	US\$ 21.64	<a href="#">+Info</a>
111429	SHOFU Ca DURA-GREEN Inverted Cone IC3 12pk - #0016	US\$ 17.31	<a href="#">+Info</a>
111468	SHOFU CA DURA-GREEN Inverted Cone IC1 12pk - #0064	US\$ 21.64	<a href="#">+Info</a>
111469	SHOFU CA DURA-GREEN Inverted Cone IC3 12pk - #0066	US\$ 21.64	<a href="#">+Info</a>
111470	SHOFU CA DURA-GREEN Inverted Cone IC5 12pk - #0068	US\$ 32.74	<a href="#">+Info</a>

### Pointed Conic Rectangle Round

SHAPE CATEGORY	KN4	KN5	KN6	KN7	PC1	PC2	RD1	RD2	RD3	RD4	RE1
SHAPE											
ISO #	080	060	100	125	030	050	030	024	036	050	062
CA (12 pk)	PN -	PN -	PN -	PN -	PN 0080	PN 0081	PN 0082	PN 0083	PN 0084	PN 0085	PN 0086
FG (12 pk)	PN -	PN -	PN -	PN -	PN -	PN 0107	PN 0108	-	-	-	-

Round

Code	Name	Price	+
111477	SHOFU CA DURA-GREEN Rectangle Re1 12pk - #0086	US\$ 16.88	<a href="#">+Info</a>
111473	SHOFU CA DURA-GREEN Round Rd1 12pk - #0082	US\$ 21.64	<a href="#">+Info</a>

111474	SHOFU CA DURA-GREEN Round Rd2 12pk - #0083	US\$ 19.26	<a href="#">+Info</a>
111475	SHOFU CA DURA-GREEN Round Rd3 12pk - #0084	US\$ 21.64	<a href="#">+Info</a>
111476	SHOFU CA DURA-GREEN Round Rd4 12pk - #0085	US\$ 21.64	<a href="#">+Info</a>
111471	SHOFU CA DURA-GREEN Point PC1 12pk - #0080	US\$ 21.64	<a href="#">+Info</a>
111472	SHOFU CA DURA-GREEN Point PC2 12pk - #0081	US\$ 19.26	<a href="#">+Info</a>

### Taper Conic

SHAPE CATEGORY	RE2	TC1	TC2	TC3	TC4	WH1	WH2	WH3	WH4	WH5	WH6	
SHAPE												
CA (12 pk)	PN	-	0089	0090	-	0092	-	-	0095	0096	-	0098
FG (12 pk)	PN	-	0109	-	-	-	0110	-	-	-	-	-

Code	Name	Price	+
111478	SHOFU CA DURA-GREEN Tapper Tc1 12pk - #0089	US\$ 17.31	<a href="#">+Info</a>
111480	SHOFU CA DURA-GREEN Tapper Tc4 12pk - #0092	US\$ 21.64	<a href="#">+Info</a>

### Wheel

Code	Name	Price	+
111481	SHOFU CA DURA-GREEN Wheel WH3 12pk - #0095	US\$ 21.64	<a href="#">+Info</a>
111483	SHOFU CA DURA-GREEN Wheel WH6 12pk - #0098	US\$ 21.64	<a href="#">+Info</a>

## Abrasive Mounted Stones

### Mounted Points Kits



### Mounted Points Kits - SHOFU

Mounted Points Kits - SHOFU

#### Ceramaster

Ceramaster

Code	Name	Price	+
111660	All Ceramic Finish/Polish Kit	US\$ 128.86	<a href="#">+Info</a>

#### Dura-Green

Dura-Green

Code	Name	Price	+
111661	All Ceramic Finishing & Polish	US\$ 156.65	<a href="#">+Info</a>

## Abrasive Mounted Stones

### RA-Mounted Stones



### RA WHITE ARKANSAS MOUNTED STONES - BESQUAL

RA WHITE ARKANSAS MOUNTED STONES - BESQUAL



Aluminum Oxide. For fine finishing of composites.

Code	Name	Price	+
102807	BESQUAL RA CN1 WHITE Mounted Stones - Cone Shape 12pk - #106-201	US\$ 9.38	+Info
102808	BESQUAL RA CY1 WHITE Mounted Stones - Cylinder Shape 12pk - #106-202	US\$ 7.90	+Info
102809	BESQUAL RA FL2 WHITE Mounted Stones - Flame Shape 12pk - #106-204	US\$ 9.57	+Info
102810	BESQUAL RA FL3 WHITE Mounted Stones - Flame Large Shape 12pk - #106-205	US\$ 7.99	+Info
102811	BESQUAL RA RD1 WHITE Mounted Stones - Round Shape 12pk - #106-207	US\$ 7.90	+Info
102812	BESQUAL RA IC1 WHITE Mounted Stones - Inverted Cone Shape 12pk - #106-206	US\$ 8.29	+Info

## RA GREEN ARKANSAS MOUNTED STONES - BESQUAL

Silicon Carbide For Porcelain, Gold, Amalgam and Enamel



Code	Name	Price	+
102801	BESQUAL RA RD1 GREEN Mounted Stones - Round Shape 12pk - Silicone carbide for po...	US\$ 7.90	+Info
102800	BESQUAL RA CN1 GREEN Mounted Stones - Cone Shape 12pk - Silicone carbide for por...	US\$ 9.87	+Info
102802	BESQUAL RA CY1 GREEN Mounted Stones - Small Cylinder Shape 12pk - Silicone carbi...	US\$ 7.99	+Info
102803	BESQUAL RA CY2 GREEN Mounted Stones - Large Cylinder Shape 12pk - Silicone carbi...	US\$ 9.87	+Info
102804	BESQUAL RA FL2 GREEN Mounted Stones - Flame Shape 12pk - Silicone carbide for po...	US\$ 9.87	+Info
102805	BESQUAL RA IC1 GREEN Mounted Stones - Invereted Cone Shape 12pk - Silicone carbi...	US\$ 7.90	+Info
102806	BESQUAL RA TC1 GREEN Mounted Stones - Taper Shape 12pk - Silicone carbide for ...	US\$ 9.87	+Info

### Abrasive Mounted Stones

## White mounted points



## WHITE MOUNTED STONES CA - SHOFU

Aluminium oxide stones with ceramic binder

Made of micrograined aluminium oxide grits, Dura-White stones can be used for very fine finishing and smoothing of various dental materials.

### Indication

### In the dental practice

Finishing and very fine working of e.g. composite, GIC (glass ionomer cements) filling materials, porcelain, enamel, e.g. smoothing crown



preparations after contouring with diamond abrasives

### In the dental laboratory

Very fine contouring and fine finishing of porcelain or composite restorations

### Range of use

Suitable for the following materials: porcelain, composite, compomer, enamel, glass ionomer cement

### Benefits:

- \* Precision running
- \* Vibration-free finishing
- \* Consistent finishing performance
- \* Economical
- \* High edge stability
- \* Long tool life
- \* Wide range of shapes

### Working end

Aluminium oxide stones (Arkansas stones) with a ceramic binder, fine grain size (micro-grain size)

### Shapes and sizes

- \* 14 shapes in total
- \* HP shank - straight handpiece: 10
- \* CA shank - contra-angle piece: 8
- \* FG shank - turbine: 9

**US DENTAL DEPOT Inc is an Authorized SHOFU Dealer in the US.**

## CA Contra Angle Stones

Code	Name	Price	+
111570	SHOFU CA Dura-WHITE Wheel Wh4 12pk - #0232	US\$ 22.30	<a href="#">+Info</a>
111569	SHOFU CA Dura-WHITE Tapper Tc1 12pk - #0230	US\$ 22.30	<a href="#">+Info</a>
111567	SHOFU CA Dura-WHITE Rectangle 12pk - #0229	US\$ 22.30	<a href="#">+Info</a>
111566	SHOFU CA Dura-WHITE Round Rd2 12pk - #0228	US\$ 22.30	<a href="#">+Info</a>
111565	SHOFU CA Dura-WHITE Round Rd 12pk - #0227	US\$ 17.84	<a href="#">+Info</a>
111564	SHOFU CA Dura-WHITE Inverted Cone IC4 12pk - #0225	US\$ 22.30	<a href="#">+Info</a>
111563	SHOFU CA Dura-WHITE Flame FI3 12pk - #0224	US\$ 22.08	<a href="#">+Info</a>
111562	SHOFU CA Dura-WHITE Flame FI2 12pk - #0223	US\$ 22.30	<a href="#">+Info</a>
111561	SHOFU CA Dura-WHITE Flame FL1 12pk - #0222	US\$ 22.30	<a href="#">+Info</a>
111560	SHOFU CA Dura-WHITE Cone Cn1 12pk - #0221	US\$ 22.08	<a href="#">+Info</a>

## Catalog Index / Indice

FG mounted stones .....	1
Green mounted points .....	3
Mounted Points Kits .....	5
RA-Mounted Stones .....	5
White mounted points .....	6