

## Stripping Tools

### Stripping Tools..



### RECIPROCATING IPR SYSTEM - SELECT DENTAL RECIPROCATING IPR SYSTEM - SELEC DENTAL



These itmes are compatible with DIGITECH Stripping tool contra angle

A safe, precise and reliable method to achieve accurate interproximal stripping. More controllable and safer for the dentist and patient than diamond discs. Provides the flexibility that diamond and carbide burs can't. Faster and more comfortable for the patient than manual hand stripping. Maximum speed 20,000 oscillations per minute.

To be used with manual handle stripping tool: 109-1603 & 109-1604.

**DO NOT work with regular contrangles**

### CEMENT EXTRACT

Code	Name	Price	+
123084	SELECT Serrated Strip for Cement Extraction each - #109-1615	US\$ 12.72	<a href="#">+Info</a>

### DOUBLE SIDED

Code	Name	Price	+
123069	SELECT Diamond Strips Double Sided Grit 15 Micron each - #109-1605	US\$ 12.72	<a href="#">+Info</a>
123070	SELECT Diamond Strips Double Sided Grit 25 Micron each - #109-1606	US\$ 12.72	<a href="#">+Info</a>
123071	SELECT Diamond Strips Double Sided Grit 40 Micron each - #109-1607	US\$ 12.72	<a href="#">+Info</a>
123072	SELECT Diamond Strips Double Sided Grit 60 Micron each - #109-1608-H	US\$ 12.72	<a href="#">+Info</a>
123073	SELECT Diamond Strips Double Sided Grit 90 Micron each - #109-1609	US\$ 12.72	<a href="#">+Info</a>

### LEFT SINGLE SIDED

Code	Name	Price	+
123074	SELECT Diamond Strips Left Single Sided Grit 15 Micron each- #109-1620	US\$ 12.72	<a href="#">+Info</a>
123075	SELECT Diamond Strips Left Single Sided Grit 25 Micron each - #109-1622	US\$ 12.72	<a href="#">+Info</a>
123076	SELECT Diamond Strips Left Single Sided Grit 40 Micron each- #109-1624	US\$ 12.72	<a href="#">+Info</a>
123077	SELECT Diamond Strips Left Single Sided Grit 60 Micron each- #109-1626	US\$ 12.72	<a href="#">+Info</a>
123078	SELECT Diamond Strips Left Single Sided Grit 90 Micron each- #109-1627	US\$ 12.72	<a href="#">+Info</a>

## RIGHT SINGLE SIDED

Code	Name	Price	+
123079	SELECT Diamond Strips Right Single Sided Grit 15 Micron each- #109-1628	US\$ 12.72	<a href="#">+Info</a>
123080	SELECT Diamond Strips Right Single Sided Grit 25 Micron each- #109-1630	US\$ 12.72	<a href="#">+Info</a>
123081	SELECT Diamond Strips Right Single Sided Grit 40 Micron each- #109-1632	US\$ 12.72	<a href="#">+Info</a>
123082	SELECT Diamond Strips Right Single Sided Grit 60 Micron each- #109-1634	US\$ 12.72	<a href="#">+Info</a>
123083	SELECT Diamond Strips Right Single Sided Grit 90 Micron each- #109-1635	US\$ 12.72	<a href="#">+Info</a>

## DENMAT STRIPPING - CERISAW

### DENMAT - CERISAW

#### Feature

Cerisaw is a tool used for interproximal finishing. It has a handle that prevents accidental slips and injury while allowing interproximal finishing with one hand.



Code	Name	Price	+
428474	DENMAT - CERI-SAW Posterior Kit: 1 Holder with 11 Blades. For Interproximal Fin...	US\$ 100.38	<a href="#">+Info</a>
428481	DENMAT - CERI-SAW Anterior Kit: 1 Holder with 11 Blades. For Interproximal Fini...	US\$ 100.38	<a href="#">+Info</a>
115173	CERI-SAW Posterior refill blades 10pk refill for Interproximal Finishing 13371	US\$ 58.33	<a href="#">+Info</a>
115174	DENMAT - CERI-SAW Anterior 10pk Blades refill for Interproximal Finishing - # 0...	US\$ 58.33	<a href="#">+Info</a>

## PR STAINLESS STEEL HANDLE - SELECT DENTAL

### PR STAINLESS STEEL HANDLE - SELECT DENTAL

100% Stainless Steel  
Knurled handles.  
To be used with the hand saw stripping tool

**Available: Straight- Angled**



Code	Name	Price	+
123067	SELECT IPR Stainless Steel Handle: Straight fits all DS & SS Saws - Each #109-1603	US\$ 39.20	<a href="#">+Info</a>
123068	SELECT IPR Stainless Steel Handle: Angled fits all DS & SS Saws - Each #109-1604	US\$ 39.20	<a href="#">+Info</a>

## Stripping Tools.. - EASY-SMILE

### Stripping Tools.. - EASY-SMILE

Code	Name	Price	+
466267	IPR IR Stripping System	US\$ 26.00	<a href="#">+Info</a>

## Stripping Tools.. - HOUSE BRAND

### Stripping Tools.. - HOUSE BRAND

Code	Name	Price	+
------	------	-------	---

122939	Mechanical Diam Stripping Syst	US\$ 135.00	<a href="#">+Info</a>
--------	--------------------------------	-------------	-----------------------

## Stripping Tools.. - MICRODONT

### Stripping Tools.. - MICRODONT

Code	Name	Price	+
470307	IPR saw blade .005 Double 10pk	US\$ 27.12	<a href="#">+Info</a>
470308	MIR 0.10mm Double Strips 10pk	US\$ 32.89	<a href="#">+Info</a>
470309	MIR 0.20mm Double Strips 10pk	US\$ 32.89	<a href="#">+Info</a>
470310	MIR 0.30mm Double strips 10pk	US\$ 32.89	<a href="#">+Info</a>
470311	MIR Double 0.40mm Strips 10pk	US\$ 32.89	<a href="#">+Info</a>
470312	IPR Double 0.50mm Purple 10pk	US\$ 32.89	<a href="#">+Info</a>
470313	MIR 0.10mm IGHT Strips 10pk	US\$ 32.89	<a href="#">+Info</a>
470314	MIR 0.20mm RIGHT Strips 10pk	US\$ 32.89	<a href="#">+Info</a>
470315	MIR 0.15mm RIGHT Strips 10pk	US\$ 32.89	<a href="#">+Info</a>
470316	MIR 0.20mm RIGHT Strips 10pk	US\$ 32.89	<a href="#">+Info</a>
470317	MIR 0.05mm LEFT Strips 10pk	US\$ 32.89	<a href="#">+Info</a>
470318	MIR 0.10mm LEFT Strips Red 10pk	US\$ 15.66	<a href="#">+Info</a>
470319	MIR 0.15mm LEFT Strips 10pk	US\$ 32.89	<a href="#">+Info</a>
470320	MIR 0.20mm LEFT Strip 10pk	US\$ 32.89	<a href="#">+Info</a>
470321	MIR 1.0 Interprox Kit 12pk + handle	US\$ 36.48	<a href="#">+Info</a>
470322	IPR Handle 10pk	US\$ 45.67	<a href="#">+Info</a>

## Stripping Tools.. - PREMIER

### Stripping Tools.. - PREMIER

Code	Name	Price	+
453281	Thin-Flex X927-7 (0.15mm, 45 micron)	US\$ 58.44	<a href="#">+Info</a>

## Stripping Tools.. - SELECT DENTAL

### Stripping Tools.. - SELECT DENTAL

Code	Name	Price	+
123047	Enamel Hand Strip Kit	US\$ 36.26	<a href="#">+Info</a>
123048	Diamond Disc Guard 7/8 22mm	US\$ 46.11	<a href="#">+Info</a>
123049	Contouring Stripping Kit	US\$ 36.02	<a href="#">+Info</a>
123050	Mini Stripper Kit	US\$ 52.92	<a href="#">+Info</a>
123051	Cello Stripper Kit	US\$ 33.08	<a href="#">+Info</a>
473696	Cello Stripper Single Sided 3mm	US\$ 27.13	<a href="#">+Info</a>

**TRIDENT STRIPPING TOOLS - TAUB**  
**TRIDENT STRIPPING TOOLS - TAUB**

These products are specializes in Oral and Surgical instruments from forceps and mouth gags to tweezers, surgical scissors and more



Code	Name	Price	+
122409	TAUB - Trident - Interproximal Saw Kit - (handle / 3x Saw / 2x Abrasive) - #99-1011	<b>US\$ 77.65</b>	<a href="#">+Info</a>
122410	TAUB - Trident - Interproximal Saw - Saw Blade (handle / Saw) - #99-1012	<b>US\$ 50.25</b>	<a href="#">+Info</a>
122411	TAUB - Trident - Interproximal Saw - Abrasive Blade (handle / Abrasive) - #99-1013	<b>US\$ 50.25</b>	<a href="#">+Info</a>
122412	TAUB - Trident - Interproximal Saw Blades Only - Replace - 5 pk - #99-1015	<b>US\$ 41.12</b>	<a href="#">+Info</a>
122413	TAUB - Trident - Interproximal Abrasive Blades Only - Replace - 5 pk - #99-1016	<b>US\$ 41.12</b>	<a href="#">+Info</a>
122414	TAUB - Trident - Interproximal Blades (3x Saw / 2x Abrasive) - 5 Pc - #99-1017	<b>US\$ 41.12</b>	<a href="#">+Info</a>

## Catalog Index / Indice

Stripping Tools..... 1