

Surgical Sutures

Sutures w/needles



NON-ABSORBABLE SUTURES

BIOVIVA Silk Suture 4-0 Reverse Cutting 19Mm 3/8" 45Cm Dozen/Box #474088#SBB683

Reliability and Quality You can Count On - Bioviva Sutures and Needles feature super fine needles and high knot comfort and security. Bioviva Sutures and Needles will enable you to satisfy the needs of your users while achieving the highest levels of performance and quality. AISI 300 Stainless Steel needles are used for improved needle integrity, sharpness, insertion smoothness, and patient comfort. Silk Black Braided Suture, ADVANTAGES:



- Super fine needles
- High knot comfort and security
- AISI 300 Stainless Steel needles
- Improved needle integrity, sharpness, insertion, smoothness and patient comfort
- 2 year tensile strength retention rating

PRODUCT DETAILS: PACKAGING 12 per Box

MONOFILAMENT (PTFE)

TFE suture is a monofilament non-absorbable, extremely soft and smooth material composed of a strand of polytetrafluoroethylene (PTFE). This material is a synthetic fluoropolymer of tetrafluoroethylene that is very strong, knots well and biologically inert.

Nylon

Nylon

HEXA

HEXA

Code	Name	Price	+
476866	NYLON Black Monofilament, Size 4/0,	US\$ 17.10	+Info

MTEC

MTEC

Code	Name	Price	+
122940	#1 Nylon Green 4-0 12pk	US\$ 13.00	+Info

SILK

3/0
3/0

BIOVIVA

BIOVIVA

Code	Name	Price	+
474086	Silk Suture 3-0 Reverse Cut-19mm	US\$ 15.40	+Info
474087	Silk Suture 3-0 Reverse Cut-22mm	US\$ 17.14	+Info
474089	BIOVIVA Silk Suture 3-0 Reverse Cut	US\$ 17.07	+Info

HEXA

HEXA

Code	Name	Price	+
476720	SILK, Size 3/0, Needle 15mm, Length	US\$ 13.78	+Info
476721	SILK, Size 3/0, Needle 19mm, Length	US\$ 13.78	+Info
476722	SILK, Size 3/0, Needle 24mm, Length	US\$ 13.78	+Info
476724	SILK, Size 3/0, Needle 24mm, Length	US\$ 13.78	+Info

MTEC

MTEC

Code	Name	Price	+
122941	#2 Silk Black 3-0 12pk	US\$ 13.00	+Info

4/0

4/0

BIOVIVA

BIOVIVA

Code	Name	Price	+
474088	BIOVIVA Silk Suture 4-0 Reverse Cut	US\$ 17.07	+Info

MTEC

MTEC

Code	Name	Price	+
122942	#3 Silk Black 4-0 12pk	US\$ 13.00	+Info

Other Models

Other Models

Code	Name	Price	+
476723	SILK, Size 4/0, Needle 19mm, Length	US\$ 13.78	+Info

5/0

5/0

Code	Name	Price	+
122944	#5 Silk Black 5-0 12pk	US\$ 12.87	+Info

REABSORVABLE SUTURES

PGA Violet Braided, Size 5/0, Needle 19mm, Length 18



Chromic Gut

The term "plain gut" is often related to a type of surgical suture material, while "PGA" stands for Polyglycolic Acid, another type of synthetic, absorbable suture material. They are not the same thing.

- [REDACTED]
- [REDACTED]

In summary, plain gut and PGA represent two different types of suture materials, with plain gut being natural and PGA being synthetic. Each has its own set of properties that make it suitable for specific types of surgical applications.

3/0
3/0

Code	Name	Price	+
474091	BIOVIVA Suture Chromic Gut C-6 3/0	US\$ 25.38	+Info
474093	BIOVIVA Suture Chromic Gut C-6 3/0	US\$ 25.38	+Info

4/0
4/0

Code	Name	Price	+
474090	BIOVIVA Suture Chromic Gut C-6 4/0	US\$ 25.38	+Info
474092	BIOVIVA Suture Chromic Gut C-6 4/0	US\$ 25.38	+Info

LIQUID GLUES

LIQUID GLUES

Code	Name	Price	+
416750	SOCKIT! Hydrogel 10gr syr	US\$ 0.00	+Info

PGA Poly Glycolic Acid

Code	Name	Price	+
------	------	-------	---

476867	HEXA - Suture PGA Violet Braided, Size 5/0, Needle 19mm, Length 18 - # HPGA-57752V	US\$ 35.10	+Info
476868	HEXA - Suture PGA Violet Braided, Size 3/0, Needle 19mm, Length 30 - # HPGA-57753U	US\$ 35.10	+Info
476869	HEXA - Suture PGA Violet Braided, Size 4/0, Needle 19mm, Length 30" - # HPGA-57754U	US\$ 35.10	+Info
476870	HEXA - Suture PGA Violet Braided, Size 4/0, Needle 19mm, Length 18" - # HPGA-57751V	US\$ 35.10	+Info
476871	HEXA - Suture PGA Violet Braided, Size 3/0, Needle 19mm, Length 18" - # HPGA-57750V	US\$ 35.10	+Info

Plain Gut

Plain Gut

3/0

3/0

Code	Name	Price	+
474095	BIOVIVA Suture Plain Gut C-6 3/0 18	US\$ 26.71	+Info
474096	BIOVIVA Suture Plain Gut C-31 3/0 1	US\$ 26.71	+Info
474097	BIOVIVA Suture Plain Gut C-7 3/0 27	US\$ 26.82	+Info

4/0

4/0

Code	Name	Price	+
474094	BIOVIVA Suture Plain Gut C-6 4/0 27	US\$ 26.71	+Info
474098	BIOVIVA Suture Plain Gut C-6 4/0 18	US\$ 26.71	+Info

Surgical Sutures

Suture Needles



Suture Needles - MTEC

Suture Needles - MTEC

Code	Name	Price	+
417305	Needle Circ Curv Cut Edg 16mm 12pk	US\$ 7.43	+Info

Catalog Index / Indice

Sutures w/needles	1
Suture Needles	4