

## X-Ray Units & Systems

### Developer



### INSTANT FILM PROCESSING CHEMISTRY GALLONS - DENTX

Excel Film Processing Chemistry

Quality chemistry ensures clearer films and an efficiently running processor.

Dental offices have trusted Dent-X's Excel Developer and Fixer for decades because it is the highest quality chemistry available.

- \* IFP Case Contains 4 One-Gallon Bottles
- \* IFP Skid Contains 18 Case
- \* Ready to Use

To rise above the rest, Excel!

Excel Chemistry Features:

- \* Superior Quality Ready-To-Use Chemistry for Film Processors
- \* Survey Rated as #1 Film Processing Developer and Fixer
- \* Environmentally Friendly, Non-Hazardous Chemistry
- \* Cleaner Working Film Processor
- \* Low Oxidation Rate Ensures a Longer Shelf Life
- \* Lower Replenishment Rate in the Processor
- \* Provides High Film Output
- \* Superb Chelating Agents for Highest Quality Films
- \* Ready-To-Use, Concentrate and Instant Film Processing Agents
- \* Bottling & Distribution Facilities in California and New York
- \* Each Case Contains 4 Gallons of Chemistry
- \* Quarterly Developer and Fixer Specials



#### IFP 18 Case

Code	Name	Price	+
114291	DENTX- Single Case IFP Developer & Fixer 18 Cases - ( Conatins - 36 Gallons of D...	US\$ 1,170.00	<a href="#">+Info</a>
114292	DENTX- Single Case IFP Developer Only - 18 Cases - ( Conatins: 72 Gallons of Dev...	US\$ 1,170.00	<a href="#">+Info</a>

#### IFP Single Case

Code	Name	Price	+
114288	DENTX- Single Case IFP Developer & Fixer - 4gal ( Conatins - 2 Gallons of Develo...	US\$ 100.49	<a href="#">+Info</a>
114289	DENTX- Single Case IFP Developer Only - ( Conatins: 4 Gallons of Developer Only	US\$ 87.00	<a href="#">+Info</a>
114290	DENTX- Single Case IFP Fixer Only - ( Conatins: 4 Gallons of Fixer Only ) #99926...	US\$ 87.00	<a href="#">+Info</a>

## Developer - AIR TECHNIQUES

Developer - AIR TECHNIQUES

Code	Name	Price	+
400653	Developer & Fixer Set 2 x 1gal.	US\$ 101.34	<a href="#">+Info</a>
400655	PeriPro Developer & Fixer Set 3x1qt	US\$ 64.73	<a href="#">+Info</a>
400658	Spray 2000 Plus 1ltBottle	US\$ 77.99	<a href="#">+Info</a>

## Developer - ECCO

Developer - ECCO

Code	Name	Price	+
110652	Chairside Darkroom	US\$ 208.00	<a href="#">+Info</a>

## Developer - FLOW

Developer - FLOW

Code	Name	Price	+
124366	Value Chem Devel & Fixer 1gal.	US\$ 41.80	<a href="#">+Info</a>

## X-Ray Units & Systems

### Film & Sensors Holders



## Film & Sensors Holders - CROSSTEX

Film & Sensors Holders - CROSSTEX

Code	Name	Price	+
400330	Sensor Slippers LG #2 yellow 50pk	US\$ 41.05	<a href="#">+Info</a>

## Film & Sensors Holders - DENTSPLY

Film & Sensors Holders - DENTSPLY

Code	Name	Price	+
469237	XCP Yellow Posterior Bite Block 25/	US\$ 37.64	<a href="#">+Info</a>
469614	XCP-ORA XPC-ORA Armonly- 1/Pk. Un	US\$ 135.19	<a href="#">+Info</a>
469560	XCP-DS Combo Kit for Suni Digital S	US\$ 0.00	<a href="#">+Info</a>
469506	XCP-DS Digital Sensor Holder Kit- C	US\$ 0.00	<a href="#">+Info</a>
469388	XCP/BAI Bite-Wing Kit for Film/PSP	US\$ 0.00	<a href="#">+Info</a>
469385	XCP Endodontic Kit refill: Aiming R	US\$ 21.50	<a href="#">+Info</a>
469376	XCP/BAI #2 Film Vertical Bite-Wing	US\$ 0.00	<a href="#">+Info</a>
469370	XCP Endodontic Kit refill: Endo Arm	US\$ 63.80	<a href="#">+Info</a>
469353	XCP/BAI #0/1 Film Horizontal Bite-W	US\$ 0.00	<a href="#">+Info</a>
469349	Snap-A-Ray Universal Film Holder- 1	US\$ 0.00	<a href="#">+Info</a>
469302	XCP Endodontic Kit- new endodontic	US\$ 72.77	<a href="#">+Info</a>
469289	XCP Evolution 2000 XCP color-coded	US\$ 0.00	<a href="#">+Info</a>
469275	XCP Blue Anterior Bite Block 25/Pk.	US\$ 45.24	<a href="#">+Info</a>

469271	SUPA Disposable Foam Film Holder- s	US\$ 36.70	<a href="#">+Info</a>
469113	XCP-DS FIT Universal Sensor Holder	US\$ 60.51	<a href="#">+Info</a>
469235	XCP/BAI #2 Film Horiz. Bite-Win 6pk	US\$ 26.25	<a href="#">+Info</a>
469227	XCP XCP/BAI Anterior Aiming Ring -	US\$ 20.79	<a href="#">+Info</a>
469223	Snap-A-Ray Universal Film Holder- 3	US\$ 0.00	<a href="#">+Info</a>
469215	XCP-DS FIT Universal Sensor Holder	US\$ 86.11	<a href="#">+Info</a>
469196	XCP XCP/BAI Posterior Aiming Ring -	US\$ 24.39	<a href="#">+Info</a>
469172	XCP XCP/BAI Anterior Arm - Blue Pro	US\$ 53.94	<a href="#">+Info</a>
469170	XCP-DS FIT Endodontic Only Kit. Dig	US\$ 115.67	<a href="#">+Info</a>
469167	XCP-DS Bitewing Aiming Ring - red-	US\$ 19.62	<a href="#">+Info</a>
469163	XCP XCP/BAI Posterior Arm - Yellow	US\$ 57.72	<a href="#">+Info</a>
469144	XCP XCP/BAI Bite Wing Arm - red Pro	US\$ 55.89	<a href="#">+Info</a>
469128	Uni-Grip 360 Universal Disposable R	US\$ 87.05	<a href="#">+Info</a>
469117	XCP-DS FIT Universal Sensor Holder	US\$ 61.06	<a href="#">+Info</a>
469114	XCP-DS FIT Universal Sensor Holder	US\$ 60.41	<a href="#">+Info</a>

### Film & Sensors Holders - DXM

Film & Sensors Holders - DXM

Code	Name	Price	+
115698	DXM Digital Sensor Possitioner	US\$ 117.47	<a href="#">+Info</a>

### Film & Sensors Holders - MARK3

Film & Sensors Holders - MARK3

Code	Name	Price	+
425593	MARK3 Cargus Bite Wing Tabs, Self-A	US\$ 8.44	<a href="#">+Info</a>

## FLOW SENSOR & FILM HOLDERS



**SENSIBLE POSITIONING SYSTEM - FLOW**

**SENSIBLE POSITIONING SYSTEM - FLOW**

**Features:**

New, improved Sensibles is a "universal" sensor positioning system that does it all. Our patent-pending ratcheting bite block locks down tight to form a custom-fit around any size sensor in either the Horizontal or Vertical position. After use, simply press the quick release button and autoclave your Sensible Bite Block.

Sensibles now feature a unique "floating" bite block.

When taking a bite wing you simply slide the bite block to the approximate mid-line of your sensor! It's quick, easy, and one bite block does it all!

Sensibles come in two sizes. Large is suited for most size 1 & 2 sensors in the vertical position. Medium is suited for most size 1 & 2 sensors in the horizontal position, and size 0 sensors in either the vertical or horizontal positions. Once you try them, you'll quickly find yourself relying on just one size Sensible for the majority of your x-rays!

Color coded sizes. Quick and easy to use. Versatile and autoclavable. Economically priced too. What a Sensible way to take x-rays!



**Kits / Set**

Kits / Set

Code	Name	Price	+
124368	Univ Sensor Pos Deluxe Starter KIT	US\$ 109.91	<a href="#">+Info</a>
124369	Univ Sensor Pos Econo Starter KIT	US\$ 49.30	<a href="#">+Info</a>

**Refills**

Refills

Code	Name	Price	+
124370	Sens Positioner Refill M 12pk	US\$ 57.27	<a href="#">+Info</a>
400112	Sens Positioner Refill L 12pk	US\$ 57.27	<a href="#">+Info</a>

### RAPID POSITIONER PSP & FILM | FLOW

#### FLOW DENTAL - RAPiD Positioner for Phosphor plate & Films

For accurate x-rays, nothing beats the RAPiD system! Our RAPiD is the original color-coded film positioner. RAPiD bite blocks feature rounded corners and are flexible for easy loading and unloading. And our tabs ensure your film or PSP plates won't fall out during x-rays. Color-coding makes assembling the correct positioner quick and easy and RAPiD is autoclavable for over 100 uses.\* RAPiD is a right angle positioning device that assures correct film/cone alignment and prevents "cone cutting". It can be used for anterior, posterior, and bite-wing x-ray exams using any brand and size (3,2,1, and 0) of dental film.



### RAPID POSITIONER

#### Anterior Kit

Anterior Kit

Code	Name	Price	+
124374	RAPiD Anterior Kit Positio 0/1	US\$ 40.60	<a href="#">+Info</a>

#### Horiz Bite Wings

Horiz Bite Wings

Code	Name	Price	+
124382	RAPiD Bite-Wing Kit w/ #2 Hori	US\$ 39.15	<a href="#">+Info</a>

#### NO-Bite Wings

NO-Bite Wings

Code	Name	Price	+
124383	RAPiD Paralleling Kit	US\$ 92.30	<a href="#">+Info</a>

#### Posterior Kit

Posterior Kit

Code	Name	Price	+
124375	RAPiD Posterior Kit Posit #2	US\$ 40.60	<a href="#">+Info</a>

#### Vertical Bite Wing

Vertical Bite Wing

Code	Name	Price	+
124381	RAPiD Bite-Wing Kit w/ #2 Vert	US\$ 40.60	<a href="#">+Info</a>

#### W/ Bite Wings

w/ Bite Wings

Code	Name	Price	+
124376	RAPiD Kit w/ Bite Wing Positio	US\$ 120.70	<a href="#">+Info</a>

#### Refill

### Anterior Bite Wing

Anterior Bite Wing

Code	Name	Price	+
400113	RAPiD Posit Anterior Bite Bloc	US\$ 34.08	<a href="#">+Info</a>

### Arms refills

Arms refills

Code	Name	Price	+
124371	Bite-Wing Positioning ARM	US\$ 27.84	<a href="#">+Info</a>
124377	RAPiD Anterior Arm 1pk	US\$ 34.80	<a href="#">+Info</a>
124378	RAPiD Posterior Arm 1pk	US\$ 34.08	<a href="#">+Info</a>
400118	RAPiD Bite Wing Arm	US\$ 27.48	<a href="#">+Info</a>

### Horiz Bite Wings

Horiz Bite Wings

Code	Name	Price	+
400115	RAPiD Bite Wing Block #2 horiz 6pk	US\$ 13.85	<a href="#">+Info</a>
400117	RAPiD Bite Wing Block #0/1 horiz 6pk	US\$ 13.85	<a href="#">+Info</a>

### Post. Bite Wing

Post. Bite Wing

Code	Name	Price	+
400114	RAPiD Posit Posterior Bite Block	US\$ 34.08	<a href="#">+Info</a>

### Rings refills

Rings refills

Code	Name	Price	+
124379	RAPiD Ant/Bite Wing Ring 1pk	US\$ 9.60	<a href="#">+Info</a>
124380	RAPiD Poster. Aiming Rings 1pk	US\$ 9.60	<a href="#">+Info</a>

### Vert Bite Wings

Vert Bite Wings

Code	Name	Price	+
400116	RAPiD Bite Wing Block #2 vertical 6pk	US\$ 29.51	<a href="#">+Info</a>

## SMART POSITIONING SYSTEM

SMART POSITIONING SYSTEM - FLOW


**Features:**  
Flow's patented SMART system makes using phosphor plates just as easy, familiar and convenient as film.  
Simply load your phosphor plate into the SMART Cassette and snap to close. The disposable cassette keeps phosphor plates safe from damage and protects them from moisture and light exposure. Place the cassette into the Universal Bite Block, attach an arm and an aiming ring, then position in the mouth.  
SMART ensures perfect plate/cone alignment for every exposure. When you're done, take the plate out, process it and discard the cassette. The SMART system, including the bite blocks, arms and rings, is completely autoclavable.



Code	Name	Price	+
124367	FLOW SMART KIT Complete starter kit (PSP) - #40854 Includes: 1 Anterior Aiming A...	<b>US\$ 102.95</b>	<a href="#">+Info</a>
124372	FLOW SMART Anterior- size #1/0 cassettes (100 PCS) (PSP) - #408500	<b>US\$ 38.68</b>	<a href="#">+Info</a>
124373	FLOW SMART Posterior- size #2 cassettes- (100 PCS) (PSP) - #408520	<b>US\$ 39.88</b>	<a href="#">+Info</a>

## PLASDENT FILM & SENSOR HOLDERS

Film & Sensors Holders - PLASDENT



## SNAP-DIGITAL

SNAP-DIGITAL

Code	Name	Price	+
423677	Snap Digital Sensor Holder 3pk	<b>US\$ 6.66</b>	<a href="#">+Info</a>
423701	Plasdent Snap Digital Sensor Holder	<b>US\$ 6.90</b>	<a href="#">+Info</a>
423860	Plasdent Snap Digital Sensor Holder	<b>US\$ 6.90</b>	<a href="#">+Info</a>
423875	Plasdent Snap Digital Sensor Holder	<b>US\$ 6.90</b>	<a href="#">+Info</a>
423878	Plasdent Snap Digital Sensor Holder	<b>US\$ 6.90</b>	<a href="#">+Info</a>

424104	Plasdent Snap Digital Sensor Holder	US\$ 6.90	<a href="#">+Info</a>
--------	-------------------------------------	-----------	-----------------------

**SNAP-FILM**  
SNAP-FILM

Code	Name	Price	+
420847	Snap A-Ray Film Holder Beige-3pk	US\$ 6.53	<a href="#">+Info</a>
423537	Plasdent Snap Film Holder BLUE 3/Pk	US\$ 6.90	<a href="#">+Info</a>
423905	Plasdent Snap Film Holder - Blue Oc	US\$ 6.90	<a href="#">+Info</a>
423981	Plasdent Snap Film Holder Light Gra	US\$ 6.83	<a href="#">+Info</a>
424058	Plasdent Snap Film Holder - Green G	US\$ 6.90	<a href="#">+Info</a>
424229	Plasdent Snap Film Holder - Black K	US\$ 6.90	<a href="#">+Info</a>

**XCP System**  
XCP System

Code	Name	Price	+
400105	Plasdent X - Ray Film Aligning Syst	US\$ 35.12	<a href="#">+Info</a>

**XIOS SENSOR HOLDER BITE- WING - SIRONA**  
XIOS SENSOR HOLDER BITE- WING - SIRONA



XIOSPLUS Intraoral Sensor holder tabs are for use with XIOSPLUS Intraoral Sensor System.

Available: Yellow- Red

Code	Name	Price	+
113781	SIRONA- XIOS Sensor Holder Bite-Wing RED 100pk Positioning (#8115465)	US\$ 93.38	<a href="#">+Info</a>
113782	SIRONA- XIOS Sensor Holder Bite-Wing YELLOW 100pk Positioning HS #8114246 SIRO...	US\$ 93.61	<a href="#">+Info</a>

**X-Ray Units & Systems**  
**Film Hangers**



**X-RAY FILM HANGER - ECCO**




Code	Name	Price	+
113704	X-RAY FILME HANGER 1 CLIP	US\$ 4.00	<a href="#">+Info</a>



113705	X-RAY FILME HANGER 4 CLIPS	US\$ 4.00	<a href="#">+Info</a>
113740	X-RAY FILME HANGER 12 CLIPS	US\$ 5.90	<a href="#">+Info</a>

### X-RAY FILM HANGERS - PLASDENT



Code	Name	Price	+
424084	Plasdent Film Carrier #0 and #00, W	US\$ 10.88	<a href="#">+Info</a>
424261	Plasdent Film Carrier #1, Black Pla	US\$ 10.88	<a href="#">+Info</a>

### X-Ray Units & Systems

## Film Mounts



### Film Mounts - DENTSPLY

Film Mounts - DENTSPLY

Code	Name	Price	+
469479	GreyLight Film Mounts Plastic/Grey	US\$ 0.00	<a href="#">+Info</a>
469490	Grip Right Film Mounts Plastic/Tran	US\$ 0.00	<a href="#">+Info</a>
469489	Grip Right Film Mounts Plastic/Tran	US\$ 65.00	<a href="#">+Info</a>
469488	Grip Right Film Mounts Plastic/Tran	US\$ 37.47	<a href="#">+Info</a>
469487	Grip Right Film Mounts Plastic/Tran	US\$ 39.00	<a href="#">+Info</a>
469486	Grip Right Film Mounts Plastic/Tran	US\$ 0.00	<a href="#">+Info</a>
469485	Grip Right Film 23625 Mounts Plasti	US\$ 92.55	<a href="#">+Info</a>
469484	GreyLight Film Mounts Plastic/Grey	US\$ 51.25	<a href="#">+Info</a>
469483	GreyLight Film Mounts Plastic/Grey	US\$ 36.32	<a href="#">+Info</a>
469482	GreyLight Film Mounts Plastic/Grey	US\$ 51.25	<a href="#">+Info</a>
469481	GreyLight Film Mounts Plastic/Grey	US\$ 86.19	<a href="#">+Info</a>
469480	GreyLight Film Mounts Plastic/Grey	US\$ 0.00	<a href="#">+Info</a>
469491	Grip Right Film Mounts Plastic/Tran	US\$ 74.19	<a href="#">+Info</a>
469478	GreyLight Film Mounts Plastic/Grey	US\$ 51.94	<a href="#">+Info</a>
469477	GreyLight Film Mounts Plastic/Grey	US\$ 51.94	<a href="#">+Info</a>
469476	GreyLight Film Mounts Plastic/Grey	US\$ 95.62	<a href="#">+Info</a>
469475	GreyLight Film Mounts Plastic/Grey	US\$ 116.73	<a href="#">+Info</a>
469474	GreyLight Film Mounts Plastic/Grey	US\$ 84.90	<a href="#">+Info</a>
469473	GreyLight Film Mounts Plastic/Grey	US\$ 89.49	<a href="#">+Info</a>
469472	GreyLight Film Mounts Plastic/Grey	US\$ 89.49	<a href="#">+Info</a>
469471	Klear-Vue Film Mounts Cardboard/Whi	US\$ 51.25	<a href="#">+Info</a>
469470	Klear-Vue Film Mounts Cardboard/Whi	US\$ 66.55	<a href="#">+Info</a>
469469	Klear-Vue Film Mounts Cardboard/Whi	US\$ 136.98	<a href="#">+Info</a>
469502	Grip Right Film Mounts Plastic/Tran	US\$ 40.12	<a href="#">+Info</a>
469584	Transview Masked translucent plasti	US\$ 38.94	<a href="#">+Info</a>
469518	EZ Tab #30-6011 gray opaque vinyl t	US\$ 0.00	<a href="#">+Info</a>
469517	EZ-View Transparent vinyl pocket wi	US\$ 0.00	<a href="#">+Info</a>
469515	EZ-View X-Ray Film Mounts- Masked P	US\$ 0.00	<a href="#">+Info</a>
469513	Eezeemount Gray cardboard die-cut s	US\$ 70.85	<a href="#">+Info</a>
469510	Klear-Vue Film Mounts Cardboard/Whi	US\$ 94.08	<a href="#">+Info</a>

469509	Klear-Vue Film Mounts Cardboard/Whi	US\$ 43.59	+Info
469508	Klear-Vue Film Mounts Cardboard/Whi	US\$ 94.08	+Info
469505	Grip Right Film Mounts Plastic/Tran	US\$ 26.39	+Info
469504	Grip Right Film Mounts Plastic/Tran	US\$ 0.00	+Info
469503	Grip Right Film Mounts Plastic/Tran	US\$ 43.59	+Info
469468	Klear-Vue Film Mounts Cardboard/Whi	US\$ 94.08	+Info
469501	Grip Right Film Mounts Plastic/Tran	US\$ 58.89	+Info
469500	Grip Right Film Mounts Plastic/Tran	US\$ 63.48	+Info
469499	Grip Right Film Mounts Plastic/Tran	US\$ 63.48	+Info
469498	Grip Right Film Mounts Plastic/Tran	US\$ 69.60	+Info
469497	Grip Right Film Mounts Plastic/Tran	US\$ 92.55	+Info
469496	Grip Right Film Mounts Plastic/Tran	US\$ 101.74	+Info
469495	Grip Right Film Mounts Plastic/Tran	US\$ 101.74	+Info
469494	Grip Right Film Mounts Plastic/Tran	US\$ 106.33	+Info
469493	Grip Right Film Mounts Plastic/Tran	US\$ 52.77	+Info
469492	Grip Right Film Mounts Plastic/Tran	US\$ 80.31	+Info
469337	Klear-Vue Film Mounts Cardboard/Whi	US\$ 44.82	+Info
469397	Klear-Vue Film Mounts Cardboard/Whi	US\$ 34.06	+Info
469395	GreyLight Film Mounts Plastic/Grey	US\$ 83.20	+Info
469392	EZ-View X-Ray Film Mounts- Masked P	US\$ 0.00	+Info
469389	EZ-View X-Ray Film Mounts- Masked P	US\$ 0.00	+Info
469381	GreyLight Film Mounts Plastic/Grey	US\$ 112.78	+Info
469375	Klear-Vue Film Mounts Cardboard/Whi	US\$ 112.88	+Info
469371	GreyLight Film Mounts Plastic/Grey	US\$ 106.33	+Info
469365	GreyLight Film Mounts Plastic/Grey	US\$ 47.76	+Info
469363	Perfect Pocket Film Mounts Vinyl/BI	US\$ 0.00	+Info
469358	EZ-View X-Ray Film Mounts- Masked P	US\$ 259.93	+Info
469339	GreyLight Film Mounts Plastic/Grey	US\$ 44.98	+Info
469402	Perfect Pocket Film Mounts Vinyl/BI	US\$ 108.91	+Info
469334	Perfect Pocket Film Mounts Vinyl/BI	US\$ 65.74	+Info
469326	Perfect Pocket Film Mounts Vinyl/BI	US\$ 52.53	+Info
469321	EZ-View X-Ray Film Mounts- Masked P	US\$ 0.00	+Info
469308	GreyLight Film Mounts Plastic/Grey	US\$ 43.20	+Info
469305	Perfect Pocket Film Mounts Vinyl/BI	US\$ 84.16	+Info
469300	Perfect Pocket Film Mounts Vinyl/BI	US\$ 107.63	+Info
469295	Perfect Pocket Film Mounts Vinyl/BI	US\$ 34.15	+Info
469290	Perfect Pocket Film Mounts Vinyl/BI	US\$ 30.95	+Info
469277	EZ-View X-Ray Film Mounts- Masked P	US\$ 0.00	+Info
469258	EZ-View X-Ray Film Mounts- Masked P	US\$ 0.00	+Info
469442	Klear-Vue Film Mounts Cardboard/Whi	US\$ 41.52	+Info
469467	Klear-Vue Film Mounts Cardboard/Whi	US\$ 94.08	+Info
469466	Klear-Vue Film Mounts Cardboard/Whi	US\$ 94.08	+Info
469465	Klear-Vue Film Mounts Cardboard/Whi	US\$ 0.00	+Info
469464	Klear-Vue Film Mounts Cardboard/Whi	US\$ 0.00	+Info
469463	Perfect Pocket Film Mounts Vinyl/BI	US\$ 0.00	+Info
469462	Perfect Pocket Film Mounts Vinyl/BI	US\$ 150.70	+Info
469461	Perfect Pocket Film Mounts Vinyl/BI	US\$ 0.00	+Info
469460	Perfect Pocket Film Mounts Vinyl/BI	US\$ 64.67	+Info
469459	Perfect Pocket Film Mounts Vinyl/BI	US\$ 136.98	+Info

469443	EZ-View X-Ray Film Mounts- Masked P	US\$ 0.00	+Info
469204	Perfect Pocket Film Mounts Vinyl/BI	US\$ 123.49	+Info
469432	Klear-Vue Film Mounts Cardboard/Whi	US\$ 66.55	+Info
469425	Klear-Vue Film Mounts Cardboard/Whi	US\$ 54.78	+Info
469422	GreyLight Film Mounts Plastic/Grey	US\$ 0.00	+Info
469419	Klear-Vue Film Mounts Cardboard/Whi	US\$ 51.92	+Info
469418	Klear-Vue Film Mounts Cardboard/Whi	US\$ 86.83	+Info
469417	Perfect Pocket Film Mounts Vinyl/BI	US\$ 0.00	+Info
469415	EZ-View X-Ray Film Mounts- Masked P	US\$ 0.00	+Info
469414	Eezeemount Gray cardboard die-cut s	US\$ 62.53	+Info
469409	GreyLight Film Mounts Plastic/Grey	US\$ 92.55	+Info
469404	Perfect Pocket Film Mounts Vinyl/BI	US\$ 50.95	+Info

### Film Mounts - FLOW

Film Mounts - FLOW

Code	Name	Price	+
469429	Klear-Vue Film Mounts Cardboard/Whi	US\$ 64.92	+Info

### Film Mounts - MARK3

Film Mounts - MARK3

Code	Name	Price	+
425456	MARK3 Cargus X-Ray Film Mounts Univ	US\$ 26.32	+Info

### Film Mounts - TPC

Film Mounts - TPC

Code	Name	Price	+
464131	TPC Advanced Technology Open Window	US\$ 11.60	+Info
464141	TPC Advanced Technology Notched Cle	US\$ 19.60	+Info
464140	TPC Advanced Technology Open Window	US\$ 0.00	+Info
464139	TPC Advanced Technology Open Window	US\$ 0.00	+Info
464138	TPC Advanced Technology Notched Cle	US\$ 16.00	+Info
464137	TPC Advanced Technology Clear pocke	US\$ 17.34	+Info
464136	TPC Advanced Technology Clear pocke	US\$ 23.20	+Info
464135	TPC Advanced Technology Open Window	US\$ 0.00	+Info
464134	TPC Advanced Technology Clear pocke	US\$ 12.80	+Info
464133	TPC Advanced Technology Notched Cle	US\$ 25.60	+Info
464132	TPC Advanced Technology Clear Pocke	US\$ 11.20	+Info
464086	TPC Advanced Technology Universal C	US\$ 38.40	+Info
464130	TPC Advanced Technology Notched Cle	US\$ 11.60	+Info
464129	TPC Advanced Technology Notched Cle	US\$ 32.00	+Info
464128	TPC Advanced Technology Open Window	US\$ 17.60	+Info
464127	TPC Advanced Technology Open Window	US\$ 20.40	+Info
464126	TPC Advanced Technology Notched Cle	US\$ 38.80	+Info
464125	TPC Advanced Technology Open Window	US\$ 0.00	+Info
464124	TPC Advanced Technology Clear pocke	US\$ 19.20	+Info

464101	TPC Advanced Technology Notched Cle	US\$ 12.51	<a href="#">+Info</a>
464100	Advanced Technology Clear pocket	US\$ 32.00	<a href="#">+Info</a>
464088	TPC Advanced Technology Clear pocke	US\$ 32.00	<a href="#">+Info</a>

## X-Ray Units & Systems

### Film Processors

## 810 ACCESSORIES - DENTX

### 810 FILM PROCESSOR ACCESSORIES DENTEX

#### Panoramic Daylight Loader :

The 810 Basic or 9000 processor may be used in a normally lighted room by employing the following optionally available Daylight Loader. This light-tight film loader accomodates all intraoral, panoramic, 5" X 7" (12.7 X 17.8 cm) and 8" X 10" (20.3 X 25.4 cm) film sizes.



#### Water Recirculator Unit:

This optional system is capable of providing wash water when it is impractical to connect the processor to permanent plumbing. It consists of a submersible pump with line cord, hose for connecting the pump to the processor and a seven gallon water container. The pump line cordplugs into the same connector as would the Solenoid Assembly.

#### Mounting Base:

The Mounting Base raises the Processor 4 1/2" (11.4 cm) above the countertop, providing clearance for all hoses and eliminating the necessity of drilling holes in the countertop.

#### Processor Stand:

The Processor Stand provides for convenient support of the Processor, Water Recirculator and storage. It is of molded construction and extremely rugged

Code	Name	Price	+
114283	DENTX- Day Light Loader #9992405001 #999-2405-001	US\$ 708.50	<a href="#">+Info</a>
114285	DENTX- 810 Plus Automatic Replenishment Upgrade #9992405029 #999-2405-029	US\$ 710.50	<a href="#">+Info</a>
114282	DENTX- Water Recirculator (includes Chemistry Drain Kit 873-020007) #9992405004 ...	US\$ 645.25	<a href="#">+Info</a>
114284	DENTX- Chemistry Drain Kit	US\$ 89.52	<a href="#">+Info</a>
114286	DENTX- Quick Silver Recovery Kit (EPA Compliance) #9992405007 #999-2405-007	US\$ 484.30	<a href="#">+Info</a>
114287	DENTX- Quick Silver Recovery Cartridge #9992405008 #999-2405-008	US\$ 456.75	<a href="#">+Info</a>

## 810 FILM PROCESSORS - DENTX

### AUTOMATIC FILM PROCESSOR

**Features:**  
There are over 65,000 Dent-X 810 automatic processors installed worldwide, spanning over 25 years of high quality film processing. Each of the three processor models share some of the same popular characteristics:

**All Dent-X Automatic Film Processors:**

- \* Are extremely low maintenance
- \* Alert the operator when the system is ready to process film
- \* Quickly and quietly produce a dry, archive-ready film
- \* Have a compact form, specifically designed for dental offices
- \* Can be cleaned and prepped in minutes
- \* Are available with an optional chemical drain kit to make cleaning even easier

### The 810 Plus and 810 Plus AR Processors

- \* Are designed for high volume offices
- \* Activate automatically once the system is ready - place the film and walk away!
- \* Shut down automatically when finished - saving energy and chemicals

\* An automatic chemical replenishment kit is included with the 810 Plus AR and optional on the 810 Plus



Code	Name	Price	+
114271	DENTX- 810 Plus Auto Film Processor ( Automatic Developer ) #9992405700 #999-240...	US\$ 4,912.70	+Info
114272	DENTX- 810 Plus AR Auto Film Processor ( Automatic Developer ) #9992405300 #999-...	US\$ 6,243.60	+Info

## X-Ray Units & Systems

### Films



## Films - CARESTREAM HEALTH

### Films - CARESTREAM HEALTH

Code	Name	Price	+
425823	Ultra-speed DF-50 #4 Occlusal film	US\$ 121.68	+Info

## INTRAORAL FILMS - KODAK EURO



Code	Name	Price	+
400745	Intraoral Film DF-57 #2 -130pk	US\$ 80.93	+Info

400746

Intraoral film DF-58 Ultrasp #2

US\$ 80.33

[+Info](#)

X-Ray Units & Systems  
**Panoramic X-Ray**



**PANOURA 18S X-RAY EQUIPMENT**

**DENTEX Panoura 18S Digital Panoramic X- Ray - 3-in-1 system**

The Panoura is a 3-in-1 panoramic system

The Panoura 18 provide yo with freedom and flexibility for in-office upgrade or 3D.

You can install the unit today and then upgrade to CEPHALOMETRIC module or with a 3D SENSOR

**CMOS SENSOR**

The Panoura 18S has a built-in cooling system located in the Sensor Housing.

The Panoura 18S also has mulitple examinations:

- \* **Panaromic 3D Imaging \***

**TMJ**

Sinus and Cephalometric

**Technicial Information:**

- \* Built-in Cooling system
- \* Multiple Examinations
- \* Autofocus Imaging
- \* Multifocal layer technology
- \* 4,500 images in a single pan
- \* CdTe technology: 3 times resolution
- \* Superior algorithm engine

For more specific technical information about PANOURA 18S system please read below .

The Panoura CMOS sensor can be easily placed on the Cephalometric arm or the panoramixc port. Just twist the lock in place and you are ready to acquire.

**EVASOFT - DENTAL IMAGING SOFTWARE**

EVAsoft Our dental imaging software supplied with the unit, is a powerful image management tool to acquire, process, print and store your images in a simple and effective way. Easy to learn and simple to operate, EVAsoft has all of the advanced features necessary to deliver rapid, effective diagnosis.

From the simplicity of tab navigation to the outstanding value it offers, EVAsoft Dental Image Managementsoftware is the choice for any practice! EVAsoft interfaces easily with all popular dental imagingsoftware suites.

**EVAsoft Features:**

- \* In-Office Training
- \* Lifetime Remote Support.
- \* 3D versions come with 3D reconstruction software.
- \* Tab Navigation



- \* Interfaces easily with all popular dental imaging software suites.
- \* Warranty: 2 Year Limited

**Panoramic X-Ray Ships  
in One Box: 92"L x 20"W x 48"H. 419 lbs. ( 200 kg )**

### PANOURA 18S PANORAMIC DIGITAL

**PANOURA 18S PANORAMIC DIGITAL Contains:**

Panoura, EVAsoft Imaging Software,  
In-Office Training,  
Lifetime Remote Support.

All Panouras may be upgraded to 3D at 2012 any time.  
3D versions come with 3D reconstruction software.

**Warranty: 2 Year Limited**



Code	Name	Price	+
114261	PANOURA 18S 2D PANORAMIC / NO CEPHALOMETIRC Upgradeable to 3D #9992433000 #999-2...	US\$ 0.00	+Info
114262	PANOURA 18S 2D PANORAMIC & CEPHALOMETRIC Upgradeable to 3D #9992433100 #999-2433...	US\$ 0.00	+Info
114263	PANOURA 18S 3D PANORAMIC #9992433300 #999-2433-300 9992433300	US\$ 0.00	+Info
114264	PANOURA 18S 3D PANORAMIC & CEPHALOMETRIC #9992433400 #999-2433-400 9992433400	US\$ 0.00	+Info

### UPGRADES FOR PANOURA 18S 2D DIGITAL

Code	Name	Price	+
114268	PANOURA 18S 2D CEPHALOMETRIC UPGRADE #9992433210 #999-2433-210 9992433210	US\$ 0.00	+Info
114269	PANOURA 18S 3D UPGRADE from 2D Panoramic #9992433313 #999-2433-313 9992433313	US\$ 0.00	+Info

### EVASOFT PANOURA SOFTWARE

Code	Name	Price	+
114266	EVAsoft Imaging Software #9992490800 - 9992490800 ( Panoura 18s Systems/ EVA Sel...	US\$ 1,300.00	+Info

### ACCESORIES

#### DENU

DENU

Code	Name	Price	+
114265	Software ProIntegrate Mgmt B	US\$ 725.00	+Info

#### IMAGEWORKS

IMAGEWORKS

Code	Name	Price	+
114267	Panoura Fr Standing Base	US\$ 2,080.00	+Info

### X-Ray Units & Systems

## Portable X-Rays



## DIOX & PROX WIRELESS PORTABLE X-RAY UNIT - DIGIMED

### DIOX & PROX WIRELESS PORTABLE X-RAY UNIT - DIGIMED

The PROX Wireless Portable X-ray is an advanced high frequency field dental x-ray Apparatus rated at 60kVDC and 2mA Tube Current which is designed to produce diagnostic quality x-ray for film and digital sensors

#### Features

- \* Optimized X-Ray Device With Digital Sensors
- \* Exposure time is adjustable by 0.01 sec. for more control
- \* Preprogrammed Buttons for fast anatomic selection with film or digital.
- \* More than 300 Times of consecutive use is available with high capacity rechargeable battery stable with P.M.C ( Protection Circuit Module ) application.
- \* Remote control available with the world's first infrared ray remote controller
- \* Since it is the only ultra slim X-Ray Device in the nation, even female users can handle without any convenience ( H.V Unit Manufacturing is Patented )
- \* Optimized Grip with non-slip silicone grip and handle.
- \* Minimized Covering since most parts including X-Ray Tube, High Voltage generator, P.C.B and Cone ( Radiation field limit device ) are shielding with lead.
- \* Can be wall mounted with optional extensions



### PROX Wireless X-Ray unit

Code	Name	Price	+
115705	DIGIMED - PROX Wireless Portable Xray - Unit - With Battery Charger 100/240volts...	US\$ 3,412.51	+Info



**X-RAYS ACCESSORIES - DIGIMED**



**Other Models**

Other Models

Code	Name	Price	+
122330	Wall Mount Arm - 88cm to 155cm	<b>US\$ 1,312.50</b>	<a href="#">+Info</a>
122331	Mobile Stand Arm 200cm - 150cm	<b>US\$ 1,367.52</b>	<a href="#">+Info</a>
122332	Portable Stand Arm	<b>US\$ 902.99</b>	<a href="#">+Info</a>
123490	BackScatter for Diox & Prox	<b>US\$ 240.00</b>	<a href="#">+Info</a>
123519	Long Cone 4 Diox	<b>US\$ 80.00</b>	<a href="#">+Info</a>
123872	Long Cone 4 Prox	<b>US\$ 80.00</b>	<a href="#">+Info</a>

**Parts for DIOX X-Ray**

**DIGIMED**

DIGIMED

Code	Name	Price	+
115768	DIOX 602 Battery DC 24V	<b>US\$ 285.00</b>	<a href="#">+Info</a>
115770	DIOX 602 Battery Charger	<b>US\$ 80.00</b>	<a href="#">+Info</a>
119566-1	DIOX Remote Control	<b>US\$ 47.04</b>	<a href="#">+Info</a>
122333	DIOX 602 pcb Control Board	<b>US\$ 352.80</b>	<a href="#">+Info</a>
420157	DIOX Wireless Remote Control	<b>US\$ 104.00</b>	<a href="#">+Info</a>

**DIOX**

DIOX

Code	Name	Price	+
124510	DIOX Power Cord 2ft 110v	<b>US\$ 9.52</b>	<a href="#">+Info</a>

**Parts for PROX X-Ray**

Code	Name	Price	+
115771	PROX Battery Charger - AC 100-240V - for Prox Portable X-Ray Unit - #115771	<b>US\$ 101.50</b>	<a href="#">+Info</a>
119566	PROX Remote Control - for Prox Portable X-Ray Unit - #	<b>US\$ 33.12</b>	<a href="#">+Info</a>
122333-1	PROX PCB Control Board Replacement for Diox 602 Portable X-Ray Unit - #	<b>US\$ 352.80</b>	<a href="#">+Info</a>

115769	PROX Battery DC 24V Replacement for Prox Portable X-Ray Unit #115769	US\$ 285.00	<a href="#">+Info</a>
--------	--	-------------	-----------------------

**Portable X-Rays - GENORAY**  
Portable X-Rays - GENORAY



**PORT-X IV**  
PORT-X IV

Code	Name	Price	+
464884	PORT-X IV Portable X-Ray	US\$ 4,597.95	<a href="#">+Info</a>
474391	PORT-X IV Portable X-Ray	US\$ 9,288.80	<a href="#">+Info</a>

**Port-XII**  
Port-XII

Code	Name	Price	+
102303	PORT-XII Portable X-Ray	US\$ 3,700.00	<a href="#">+Info</a>
123530	ZEN-PX2 Portable X-Ray	US\$ 3,700.00	<a href="#">+Info</a>

**X-RAYS ACCESSORIES - GENORAY**



Code	Name	Price	+
475625	Long Cone ONLY Genoray	US\$ 0.01	<a href="#">+Info</a>

**Portable X-Rays - MAXRAY**

Portable X-Rays - MAXRAY



Code	Name	Price	+
123512	MAX-RAY Portable X-ray Blk	US\$ 3,967.50	<a href="#">+Info</a>
475200	MAX-RAY Cocoon X-RAY WHITE	US\$ 5,625.00	<a href="#">+Info</a>
123513	MAX-RAY Portable X-ray Wht	US\$ 3,967.50	<a href="#">+Info</a>
470705	MAX-RAY 5 years Warranty	US\$ 783.39	<a href="#">+Info</a>

### Portable X-Rays - WOODPECKER

Portable X-Rays - WOODPECKER

Code	Name	Price	+
474420	AI RAY PRO - Portable X-Ray Unit	US\$ 2,812.00	+Info
474915	SCATTER SHIELD for - AI RAY PRO	US\$ 292.50	+Info

### Portable X-Rays - XPECT VISION

Portable X-Rays - XPECT VISION

Code	Name	Price	+
476897	Xpect Vision Sensor #1	US\$ 2,282.00	+Info
476898	Xpect Vision Sensor #1	US\$ 4,491.90	+Info

### X-Ray Units & Systems

## Protectors /Aprons /Ponchos



### LEAD X-RAY PROTECTIVE APRONS - FLOW DENTAL

#### FLOW - LEAD X-RAY PROTECTIVE APRONS

#### Safe'n'Sure Protective Aprons

Flow's lead aprons provide traditional protection for your patients in a variety of styles. They're available in 7 fashionable colors to match your decor and create a relaxing environment for your patients and yourself. Considered the industry standard, Flow aprons are manufactured from lead rubber with an attractive vinyl front and a tear-resistant cling back. They clean easily with soap and water and withstand years of heavy use.

Made for years of heavy use, **all lead aprons provide 0.3mm lead equivalent protection.**

**Aprons are available with and without attached thyroid collars**, and thyroid collars can also be ordered separately. Select from any of the fabric colors shown below.

**Size: Adult - 24" x 26"**



#### Apron with Attached Collar

Code	Name	Price	+
124421	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in PURPLE - #750...	US\$ 165.92	+Info
124422	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in TEAL - #75051-14	US\$ 165.92	+Info

124423	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in MEDIUM BLUE - ...	US\$ 165.92	+Info
124424	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in BURGUNDY - #7...	US\$ 165.92	+Info
124425	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in GREY - #75051-19	US\$ 176.80	+Info
124426	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in ROYAL BLUE - ...	US\$ 165.92	+Info
124427	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in NAVY BLUE - #...	US\$ 165.92	+Info
124428	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in BEIGE - #7505...	US\$ 165.92	+Info
124429	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in BLACK - #7505...	US\$ 165.92	+Info
124430	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in TYE DYE - #75...	US\$ 165.92	+Info
124431	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in DOG PAWS - #7...	US\$ 165.92	+Info
124432	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in RED - #75051-35	US\$ 165.92	+Info
124433	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in HUNTER GREEN ...	US\$ 165.92	+Info
124434	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in SHIMMERING BR...	US\$ 165.92	+Info
124435	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in SHIMMERING GR...	US\$ 165.92	+Info
124436	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in SHIMMERING EL...	US\$ 165.92	+Info
124437	FLOW LEAD PROTECTIVE ADULT BIB APRON w/ COLLAR .3mm LEAD 24inx26in ANTIMICROBIAL...	US\$ 165.92	+Info

### Apron without Collar

Code	Name	Price	+
124387	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in PURPLE - #75050-12	US\$ 136.00	+Info
124388	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in TEAL - #75050-14	US\$ 136.00	+Info
124389	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in MEDIUM BLUE - #75050-16	US\$ 136.00	+Info
124390	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in BURGUNDY - #75050-18	US\$ 136.00	+Info
124391	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in GREY - #75050-19	US\$ 136.00	+Info
124392	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in ROYAL BLUE - #75050-20	US\$ 136.00	+Info
124393	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in NAVY BLUE - #75050-22	US\$ 136.00	+Info
124394	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in BEIGE - #75050-23	US\$ 136.00	+Info
124395	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in BLACK - #75050-25	US\$ 145.52	+Info
124396	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in TYE DYE - #75050-30	US\$ 136.00	+Info
124397	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in DOG PAWS - #75050-34	US\$ 136.00	+Info
124398	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in RED - #75050-35	US\$ 136.00	+Info
124399	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in HUNTER GREEN - #75050-36	US\$ 136.00	+Info

124400	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in SHIMMERING BRONZE - #75...	US\$ 136.00	+Info
124401	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in SHIMMERING GREEN - #750...	US\$ 136.00	+Info
124402	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in SHIMMERING ELECTRIC BLU...	US\$ 136.00	+Info
124403	FLOW LEAD PROTECTIVE ADULT BIB APRON .3mm LEAD 24inx26in ANTIMICROBIAL NAVY BLUE...	US\$ 136.00	+Info

## NO-LEAD PROTECTIVE APRONS - FLOW DENTAL

### FLOW - NO-LEAD X-RAY PROTECTIVE APRONS



Our new aprons are made of environmentally-friendly, lead-free, non-toxic materials.

They're lighter and more flexible than conventional lead aprons, but have the same protective qualities! All Flow **No-Lead Protective Aprons provide 0.25mm lead equivalency protection.** Aprons are backed with a tear-resistant cling fabric.

Aprons are available with and without attached thyroid collars, and thyroid collars can also be ordered separately.

**Size: Adult - 24" x 26"**

### Apron without Collar

Code	Name	Price	+
124404	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in PURPL...	US\$ 231.20	+Info
124405	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in TEAL ...	US\$ 231.20	+Info
124406	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in MEDIU...	US\$ 231.20	+Info
124407	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in BURGU...	US\$ 231.20	+Info
124408	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in GREY ...	US\$ 231.20	+Info
124409	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in ROYAL...	US\$ 231.20	+Info
124410	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in NAVY ...	US\$ 231.20	+Info
124411	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in BEIGE...	US\$ 231.20	+Info
124412	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in BLACK...	US\$ 231.20	+Info
124413	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in TYE D...	US\$ 231.20	+Info
124414	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in DOG P...	US\$ 231.20	+Info
124415	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in RED -...	US\$ 231.20	+Info

124416	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in HUNTE...	US\$ 231.20	+Info
124417	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in SHIMM...	US\$ 231.20	+Info
124418	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in SHIMM...	US\$ 231.20	+Info
124419	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in SHIMM...	US\$ 231.20	+Info
124420	FLOW LEAD FREE PROTECTIVE BIB APRON WITHOUT COLLAR .25mm NO-LEAD 24inx26in ANTIM...	US\$ 231.20	+Info

### Aprons with Attached Collar

Code	Name	Price	+
124438	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124439	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124440	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124441	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124442	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124443	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124444	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124445	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124446	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124447	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124448	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124449	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124450	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124451	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124452	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124453	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info
124454	FLOW LEAD FREE PROTECTIVE BIB APRON WITH ATTACHED COLLAR .25mm NO-LEAD 24inx26in...	US\$ 311.21	+Info

### Protectors /Aprons /Ponchos - BUFFALO

Protectors /Aprons /Ponchos - BUFFALO

Code	Name	Price	+
424604	Soft Disposable Blue Large 26in	US\$ 41.17	+Info

## Protectors /Aprons /Ponchos - FLOW

Protectors /Aprons /Ponchos - FLOW

Code	Name	Price	+
124353	Apron Uni-Hangar	US\$ 41.33	<a href="#">+Info</a>
419977	Childs Bib Aprons w/collar Blue	US\$ 136.32	<a href="#">+Info</a>
421063	E Wrap Thyroid Collar No Lead Gray	US\$ 457.50	<a href="#">+Info</a>

## Protectors /Aprons /Ponchos - HOUSE BRAND

Protectors /Aprons /Ponchos - HOUSE BRAND

Code	Name	Price	+
470640	Xray Apron w/ Collar Grey	US\$ 196.35	<a href="#">+Info</a>

## THYROID COLLAR LEAD-FREE - FLOW DENTAL

### FLOW - THYROID COLLAR LEAD-FREE

Our lead-substitute material provides 0.25mm lead equivalency protection while making our aprons 25% lighter than traditional aprons. And No-Lead aprons don't require special disposal.

They are better for you, better for your patients, and better for the environment! The safety of lead protection without the weight... or the lead!



**Size:** 5" x 27"

**Note:** Picture by way of illustration

Code	Name	Price	+
124455	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in PURPLE - #75056N...	US\$ 52.92	<a href="#">+Info</a>
124456	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in TEAL - #75056NL-14	US\$ 52.92	<a href="#">+Info</a>
124457	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in MEDIUM BLUE - #7...	US\$ 52.92	<a href="#">+Info</a>
124458	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in BURGUNDY - #7505...	US\$ 52.92	<a href="#">+Info</a>
124459	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in GREY - #75056NL-19	US\$ 52.92	<a href="#">+Info</a>
124460	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in ROYAL BLUE - #75...	US\$ 52.92	<a href="#">+Info</a>
124461	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in NAVY BLUE - #750...	US\$ 52.92	<a href="#">+Info</a>
124462	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in BEIGE - #75056NL-23	US\$ 52.92	<a href="#">+Info</a>
124463	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in BLACK - #75056NL-25	US\$ 55.82	<a href="#">+Info</a>
124464	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in TYE DYE - #75056...	US\$ 52.92	<a href="#">+Info</a>
124465	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in DOG PAWS - #7505...	US\$ 52.92	<a href="#">+Info</a>



124466	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in RED - #75056NL-35	US\$ 52.92	<a href="#">+Info</a>
124467	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in HUNTER GREEN - #...	US\$ 52.92	<a href="#">+Info</a>
124468	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in SHIMMERING BRONZ...	US\$ 52.92	<a href="#">+Info</a>
124469	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in SHIMMERING GREEN...	US\$ 52.92	<a href="#">+Info</a>
124470	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in SHIMMERING ELECT...	US\$ 52.92	<a href="#">+Info</a>
124471	FLOW LEAD FREE PROTECTIVE THYROID COLLAR .25mm NO-LEAD 5inx27in ANTIMICROBIAL NA...	US\$ 52.92	<a href="#">+Info</a>

### X-Ray Units & Systems

## X-Ray Accessories



## RADIATION MONITORING



## ICCARE MONITORING

### BADGE X-RAY MONITORING SERVICE - ICCARE

BADGE X-ray Monitoring Service - ICCARE

**Features:**  
The I.C. Care Radiation Exposure Monitoring Service is a NVLAP approved dosimetry service that provides early warning of X-ray problems, staff protection from accumulated exposure, high standards of image exposure and avoidance of noncompliance fines and penalties



Code	Name	Price	+
124317	ICCARE Radiation Badges Student 6 month X-ray monitoring service-1 Radiation Bad...	US\$ 58.90	<a href="#">+Info</a>
124315	ICCARE Radiation Badges Quarterly X-ray monitoring service 4 Radiation Badges s...	US\$ 137.46	<a href="#">+Info</a>
124316	ICCARE Radiation Badges Fetal 7 month X-ray monitoring service- 7 Radiation Badg...	US\$ 152.47	<a href="#">+Info</a>
124309	ICCARE Radiation Badges monthly X-ray monitoring Service x 1 year service ( Radia...	US\$ 208.44	<a href="#">+Info</a>

### INSTADOSE USB RADIATION MONITORS - ICCARE

INSTADOSE USB - ICCARE

**Features:**

Protect yourself and your staff from the possibility of excessive radiation dosage with this simple to use dosimeter badge. No more waiting to know your exposure!

Instadose Badge X-ray Monitoring Service Plug & Read USB X-Ray Monitoring Technology. Wear Instadose Badge, plug into your computer, visit the online site to record your dose and see your report without wasting any paper. Requires USB Port, Internet Service and 24MB space for software download. Monitoring service meets the requirements of every state



### Other Models

#### Other Models

Code	Name	Price	+
124308	INSTADOSE USB Green 1 year Svc	US\$ 183.87	+Info
124311	INSTADOSE USB Blue 1 year Svc	US\$ 183.87	+Info
124312	INSTADOSE USB Pink 1 year Svc	US\$ 183.87	+Info
124313	INSTADOSE USB Silver1 year Svc	US\$ 183.87	+Info
124314	INSTADOSE USB Black 1 year Svc	US\$ 183.87	+Info

### YEAR ACTIVATION



Code	Name	Price	+
475006	ICCARE - INSTADOSE USB YEAR RENEWAL monitoring service # RENEWAL.IN	US\$ 183.87	+Info

### RING RADIATION MONITORS - ICCARE

FINGER RING RADIATION MONITORS - ICCARE

**Features:**

Ring Dosimeter  
 Measures exposures to the extremities and skin from external sources of gamma radiation and x-ray exposure or beta radiation with maximum energy greater than 0.7 MeV with a minimum reportable dose of 20 mrem.  
 \* Required when working with more than 0.5 mCi of P-32, or 1 mCi of I-125, Cr  
 \* 51 or x-ray diffraction unit.  
 \* Design to wear inside of disposable gloves

\* Comes in Small, Medium, Large, and XL



### 12 Month Service

Code	Name	Price	+
124310	ICCARE Right Hand Small Finger Ring 12 Month Service MRSR	US\$ 268.78	+Info
124318	ICCARE Left Hand Small Finger Ring 12 Monthly Service MLSR	US\$ 268.78	+Info

124319	ICCARE Right Hand Medium Finger Ring 12 Monthly Service MRMR	US\$ 268.78	<a href="#">+Info</a>
124320	ICCARE Left Hand Medium Finger Ring 12 Monthly Service MLMR	US\$ 268.78	<a href="#">+Info</a>
124321	ICCARE Right Hand Large Finger Ring 12 Monthly Service MRLR	US\$ 268.78	<a href="#">+Info</a>
124322	ICCARE Left Hand Large Finger Ring 12 Monthly Service MLLR	US\$ 268.78	<a href="#">+Info</a>

### 4 Month Service

Code	Name	Price	+
124323	ICCARE Right Hand Small Finger Ring Quarterly service 4/ry QRSR	US\$ 192.24	<a href="#">+Info</a>
124324	ICCARE Left Hand Small Finger Ring Quarterly service 4/ry QLSR	US\$ 192.24	<a href="#">+Info</a>
124325	ICCARE Right Hand Medium Finger Ring Quarterly service 4/ry QRMR	US\$ 192.24	<a href="#">+Info</a>
124326	ICCARE Left Hand Medium Finger Ring Quarterly service 4/ry QLMR	US\$ 192.24	<a href="#">+Info</a>
124327	ICCARE Right Left Hand Large Finger Ring Quarterly service 4/ry QRLR	US\$ 192.24	<a href="#">+Info</a>
124328	Left Hand Large Finger Ring Quarterly service 4/ry QLLR	US\$ 192.24	<a href="#">+Info</a>

### WOLF MONITORING

#### X-Ray Accesories - WOLF MEDICAL

X-Ray Accesories - WOLF MEDICAL

Code	Name	Price	+
474413	RADIATION CONTROL Badge	US\$ 170.00	<a href="#">+Info</a>

### X-Ray Accesories - AIR TECHNIQUES

X-Ray Accesories - AIR TECHNIQUES

Code	Name	Price	+
400654	Formula 2000 Plus Cleaner 2pk	US\$ 73.37	<a href="#">+Info</a>

### X-Ray Accesories - ECCO

X-Ray Accesories - ECCO

Code	Name	Price	+
420841	Moblle Lead Porta Shield 24inx 48in	US\$ 729.00	<a href="#">+Info</a>

### X-RAYS ACCESSORIES - DENTX



Code	Name	Price	+
114276	Endos ACP Remote PushButton	US\$ 246.50	<a href="#">+Info</a>
114277	Endos ACP Remote Timer Kit	US\$ 314.65	<a href="#">+Info</a>

### X-Ray Units & Systems

## X-Ray Bite Blocks



### DISPOSABLE X-RAYS BITE BLOCKS - MARK 3



Code	Name	Price	+
425661	MARK3 Cargus Disposable X-Ray Bite	US\$ 37.73	<a href="#">+Info</a>
425554	MARK3 Cargus Disposable Bite Blocks	US\$ 37.73	<a href="#">+Info</a>
425602	MARK3 Cargus Disposable X-Ray Bite	US\$ 37.73	<a href="#">+Info</a>
425651	MARK3 Cargus Disposable X-Ray Bite	US\$ 26.93	<a href="#">+Info</a>
425660	MARK3 Cargus Disposable X-ray Bite	US\$ 26.93	<a href="#">+Info</a>
425662	MARK3 Cargus Disposable X-Ray Bite	US\$ 26.93	<a href="#">+Info</a>
425785	MARK3 Cargus Disposable X-Ray Bite	US\$ 0.00	<a href="#">+Info</a>

### DISPOSABLE X-RAYS BITE BLOCKS - PLASDENT



Code	Name	Price	+
423773	Block-4-Bite #7020 Disposable Bite	US\$ 21.09	<a href="#">+Info</a>
423751	Block-4-Bite #7040 Disposable Bite	US\$ 21.30	<a href="#">+Info</a>
423672	Block-4-Bite #7060 Disposable Bite	US\$ 21.30	<a href="#">+Info</a>
423491	Block-4-Bite Disposable X-Ray Bite	US\$ 21.30	<a href="#">+Info</a>
423576	Block-4-Bite Disposable X-Ray Bite	US\$ 21.30	<a href="#">+Info</a>
423697	Block-4-Bite Disposable X-Ray Bite	US\$ 21.30	<a href="#">+Info</a>
423765	Block-4-Bite Disposable X-Ray Bite	US\$ 21.30	<a href="#">+Info</a>
423853	Block-4-Bite Disposable X-Ray Bite	US\$ 21.30	<a href="#">+Info</a>
423866	Block-4-Bite Disposable X-Ray Bite	US\$ 21.30	<a href="#">+Info</a>
423954	Block-4-Bite Disposable X-Ray Bite	US\$ 21.30	<a href="#">+Info</a>
424049	Block-4-Bite Disposable X-Ray Bite	US\$ 21.30	<a href="#">+Info</a>
424262	Block-4-Bite Disposable X-Ray Bite	US\$ 21.30	<a href="#">+Info</a>

### X-Ray Bite Blocks - DENTSPLY

X-Ray Bite Blocks - DENTSPLY

Code	Name	Price	+
469200	Stabe Disposable Bite Blocks 100/Bo	US\$ 0.00	<a href="#">+Info</a>

**X-Ray Units & Systems**  
**X-Ray Plates Scanners**



**Phosphopr Plates Scanner Units**



**APIXIA**

APIXIA

Code	Name	Price	+
473592	SPS Scanner System	US\$ 5,499.90	<a href="#">+Info</a>

**TYRIS**

TYRIS

Code	Name	Price	+
417385	TR-100 Phosphorous Plate scanner	US\$ 3,622.40	<a href="#">+Info</a>

**WOODPECKER**

WOODPECKER

Code	Name	Price	+
474436	Free Scan - image plate scanner	US\$ 3,000.00	<a href="#">+Info</a>
474437	I SCAN - image plate scanner	US\$ 3,499.92	<a href="#">+Info</a>

**X-Ray Plates Scanners - TPC**

X-Ray Plates Scanners - TPC

Code	Name	Price	+
469713	X-RAY FILM READER (AV output)	US\$ 240.00	<a href="#">+Info</a>

469714	X-RAY FILM READER (USB)	US\$ 240.00	<a href="#">+Info</a>
--------	-------------------------	-------------	-----------------------

## X-Ray Phosphor Plates



### Air Techniques

Air Techniques

Code	Name	Price	+
400657	ScanX Intraoral Phosphor Plates #0	US\$ 106.93	<a href="#">+Info</a>

### APIXIA

APIXIA

Code	Name	Price	+
473593	PSP Size #0 4pk	US\$ 125.00	<a href="#">+Info</a>
473594	PSP Size #2 4pk	US\$ 125.00	<a href="#">+Info</a>
473595	PSP Size #1 4pk	US\$ 125.00	<a href="#">+Info</a>
473596	PSP Size #3 4pk	US\$ 152.00	<a href="#">+Info</a>

### MARK3

MARK3

Code	Name	Price	+
425511	Intraoral Phosphor Plate #2	US\$ 133.24	<a href="#">+Info</a>
425555	Intraoral Phosphor Plate #0	US\$ 70.88	<a href="#">+Info</a>
425708	Intraoral Phosphor Plate #1	US\$ 70.88	<a href="#">+Info</a>
425709	Intraoral Phosphor Plate #4	US\$ 115.02	<a href="#">+Info</a>

### Phosphor Plates - WOODPECKER

Code	Name	Price	+
474477	Phosphor Plate #0	US\$ 40.00	<a href="#">+Info</a>
474478	Phosphor Plate #1	US\$ 40.00	<a href="#">+Info</a>
474479	Phosphor Plate #2	US\$ 40.00	<a href="#">+Info</a>
474480	Phosphor Plate #3	US\$ 40.00	<a href="#">+Info</a>

### SCANX INTRAORAL PHOSHOR PLATES - AIR TECHNIQUE

S

#### SCANX INTRAORAL PHOSHOR PLATES - AIR TECHNIQUES

PSPs can comfortably accommodate children, adults, those with TMJ

issues, large tori, gaggers and elderly patients

ScanX PSPs are wireless, flexible and 30 times thinner than direct wired sensors}

The entire PSP surface area is active which means you capture nearly

40% more anatomy than the comparable sized wired sensor

ScanX PSPs do not require annual insurance

Operate in normal room lighting



**Available Size: #0- #2**

Code	Name	Price	+
420175	Air Techniques - ScanX Intraoral Phosphor Plates Size #1 - 2/pk - Air Techniques ...	US\$ 127.97	+Info
400656	Air Techniques - ScanX Intraoral Phosphor Plates Size #2 - 4/pk - Air Techniques ...	US\$ 311.64	+Info

### Visionary Image Phosphor Plates AIR TECHNIQUES - 3D De nt



Code	Name	Price	+
473979	VISIONARY Phosphor Imaging Plates C	US\$ 67.50	+Info
473980	VISIONARY Phosphor Imaging Plates C	US\$ 67.50	+Info
473981	VISIONARY Phosphor Imaging Plates C	US\$ 126.00	+Info
473982	VISIONARY Phosphor Imaging Plates C	US\$ 117.00	+Info
473983	VISIONARY Phosphor Imag. Plates #4	US\$ 160.42	+Info

### Visionary Image Phosphor Plates GENDEX - 3D Dent

Code	Name	Price	+
474062	VISIONARY Phosphor Imaging Plates C	US\$ 69.00	+Info
474063	VISIONARY Phosphor Imaging Plates C	US\$ 69.00	+Info
474064	VISIONARY Phosphor Imaging Plates C	US\$ 129.00	+Info
474065	VISIONARY Phosphor Imaging Plates C	US\$ 117.00	+Info

474066	VISIONARY Phosphor Imaging Plates C	US\$ 172.50	<a href="#">+Info</a>
--------	-------------------------------------	-------------	-----------------------

### Visionary Image Phosphor Plates SOREDEX - 3D Dent

Code	Name	Price	+
474099	Phosphor Imaging Plates S #0	US\$ 69.00	<a href="#">+Info</a>
474100	VISIONARY Phosphor Imaging Plates C	US\$ 69.00	<a href="#">+Info</a>
474101	Phosphor Imag Plates Type #2-4pk	US\$ 147.00	<a href="#">+Info</a>

## X-Ray Units & Systems

### X-Ray Sensors



## EVA DIGITAL X-RAY SYSTEMS - DENTX

### EVA DIGITAL X-RAY SYSTEMS - DENTX

Backed by the proven reliability of tens of thousands of installations, EVA provides more than just great value. This classically styled, rugged sensor is just a few millimeters thin, but delivers crisp, clear radiographs using virtually any intraoral X-ray unit.

EVA's value makes it an ideal choice for any practice - from large installations to even the most value-conscious practices. No other sensor delivers EVA's combination of superior image quality, ease of use and warranty assurance.

- \* Durable CMOS technology provides superior images using only a fraction of the power required by other digital sensor technologies
- \* Impact resistant housings are designed to minimize damage from a biting patient
- \* Sensor size #1 and #2 are similar to traditional films
- \* Hot swapping sensors between multiple computers offers greater speed, convenience and value
- \* Service multiple operatories with one sensor using the EVA docking station
- \* Small, inexpensive docking stations fit in the palm of your hand and make expansion easy
- \* 2-Year warranty for peace of mind

**Made with pride in the USA**



### Complete System

Code	Name	Price	+
114227	IMAGEWORKS - Complete System EVA Select #1 System Kits #9992412100 #999-2412-100...	US\$ 0.00	<a href="#">+Info</a>
114228	IMAGEWORKS - Complete System EVA Select #2 / Kit #9992412200 #999-2412-200 99924...	US\$ 0.00	<a href="#">+Info</a>



### Sensor Only

Code	Name	Price	+
114229	EVA Select #1 Sensor Only #9992412105 #999-2412-105 9992412105 ( Panoura 18s Sys...	US\$ 0.00	<a href="#">+Info</a>
114230	EVA Select #2 Sensor Only #9992412205 #999-2412-205 9992412205 ( Panoura 18s Sys...	US\$ 0.00	<a href="#">+Info</a>

### Software



Code	Name	Price	+
114232	ProlImage Database Porting Service (note: must be scheduled) #9992810105 #999-281...	US\$ 1,950.00	<a href="#">+Info</a>

### Accessories

Code	Name	Price	+
114233	Docking Station Kit EVA Select #9992412040 #999-2412-040 9992412040 ( Panoura 18s...	US\$ 58.00	<a href="#">+Info</a>
114234	EVA Non-Sterile Sheaths9 500/Box) 870-000256 ( Panoura 18s Systems/ EVA Select)	US\$ 326.25	<a href="#">+Info</a>
114236	RINN Positioning Aids for EVA Select 870-000310 ( Panoura 18s Systems)	US\$ 217.50	<a href="#">+Info</a>

### Loyalty Program

Code	Name	Price	+
114237	PANOURA 18S 2D PANORAMIC Upgradeable to Cephalometric and/ or 3D #9992433200 #99...	US\$ 0.00	<a href="#">+Info</a>

## X-Ray Sensor - DENT-X

X-Ray Sensor - DENT-X

Code	Name	Price	+
428527	EVA Dent-X Digital X-Ray Sensor - F	US\$ 0.00	<a href="#">+Info</a>

## X-RAY SENSOR - DENTERPRISE



### INSTA RAY

Code	Name	Price	+
420338	InstaRay USB Sensor Size 1	US\$ 3,186.15	<a href="#">+Info</a>

420339	InstaRay USB Sensor Size 2	US\$ 4,357.43	<a href="#">+Info</a>
--------	----------------------------	---------------	-----------------------

### QUICKRAY

Code	Name	Price	+
420336	QuickRay X-Ray Sensor #1	US\$ 3,622.58	<a href="#">+Info</a>
420337	QuickRay X-Ray Sensor #2	US\$ 4,194.61	<a href="#">+Info</a>
473790	QuickRay Twain Driver	US\$ 183.75	<a href="#">+Info</a>

### X-Ray Sensor - DIGIMED

X-Ray Sensor - DIGIMED

Code	Name	Price	+
115696	DVS-200 Intraoral Sensor 1.5	US\$ 3,673.00	<a href="#">+Info</a>

### X-Ray Sensor - GENORAY

X-Ray Sensor - GENORAY



### Other Models

Other Models

Code	Name	Price	+
433968	PortView Digital Sensor #1 w/ Softw	US\$ 3,900.12	<a href="#">+Info</a>
466320	PortView Digital Sensor #2 w/ Softw	US\$ 3,900.12	<a href="#">+Info</a>

### Port-View #2

Port-View #2

Code	Name	Price	+
474920	PortView Digital Sensor #2	US\$ 4,495.50	<a href="#">+Info</a>

### X-Ray Sensor - GLOBALTEC

X-Ray Sensor - GLOBALTEC

Code	Name	Price	+
421119	GLOBALTEC HDR 1000 Sensor	US\$ 2,362.50	<a href="#">+Info</a>

### X-Ray Sensor - HANDY

X-Ray Sensor - HANDY

Code	Name	Price	+
416695	HANDY HDR-500A Sensor	US\$ 1,899.98	<a href="#">+Info</a>
466389	HANDY HDR-600 Sensor	US\$ 2,362.50	<a href="#">+Info</a>

### X-Ray Sensor - IMAGEWORKS

X-Ray Sensor - IMAGEWORKS

Code	Name	Price	+
470913	EVA Select #1 PEDO Sensor Only	US\$ 0.00	<a href="#">+Info</a>

### X-Ray Sensor - WOODPECKER

X-Ray Sensor - WOODPECKER

Code	Name	Price	+
474418	H1 - Intraoral Sensor	US\$ 1,826.00	<a href="#">+Info</a>
474419	H2 - Intraoral Sensor	US\$ 2,373.80	<a href="#">+Info</a>

### X-Ray Sensor - XPECT VISION

X-Ray Sensor - XPECT VISION

Code	Name	Price	+
474407	Xpect Vision Sensor #0.8	US\$ 2,250.00	<a href="#">+Info</a>
474408	Xpect Vision Sensor #1.5	US\$ 2,249.99	<a href="#">+Info</a>

## X-Ray Units & Systems

### X-Ray Units



## ENDO ACP INTRAORAL XRAY IMAGEWORKS - DENTX

### ENDO ACP INTRAORAL X-RAY IMAGEWORKS - DENTX

#### Features:

ENDO ACP Intraoral X-Ray Unit  
Endos ACP: Optimized Parameters

The sophisticated Endos ACP timer allows the selection of optimum exposure parameters using pre-programmed anatomical techniques.

The choice of the diagnostic target is easily achieved in two quick steps: Choose the size of the patient from the three available and choose the dentition area to be examined. At the touch of a key a specific modality may be chosen which optimizes exposure times for use with any digital acquisition system.

#### Advantage:

- \* DLC Technology - (Dual Low Current)
- \* Customizable - ( Endos can be configured with a remote timer )
- \* Customizable - ( Endos tailor-fits every operator )
- \* Safety - (advanced selfdiagnostics which alert the operator to any problems that may arise )
- \* Hygiene - ( glossy finish for quick and easy cleaning to maintain the hygiene of your working environment )
- \* New Arm - (specifically designed for easy positioning )



**DENTX**

DENTX

Code	Name	Price	+
114279	Stand Only for X-Ray Endo Mobile	<b>US\$ 840.44</b>	<a href="#">+Info</a>

**IMAGEWORKS**

IMAGEWORKS

Code	Name	Price	+
417351	IMAGE SCAN HD Intraoral X- Ray	<b>US\$ 3,445.00</b>	<a href="#">+Info</a>
474311	IMAGE SCAN HD Intraoral X- Ray	<b>US\$ 3,445.00</b>	<a href="#">+Info</a>
474312	IMAGE SCAN HD Intraoral X- Ray	<b>US\$ 3,445.00</b>	<a href="#">+Info</a>

**X-Ray Units - EVOLUTION**

X-Ray Units - EVOLUTION

Code	Name	Price	+
400000	EVOLUTION X-ray Stand 110V	<b>US\$ 3,039.77</b>	<a href="#">+Info</a>
400001	EVOLUTION X-ray Wall Mnt 110V	<b>US\$ 2,414.99</b>	<a href="#">+Info</a>

**X-Ray Units & Systems**

**X-ray viewer**



**X-ray viewer - FLOW**

X-ray viewer - FLOW

Code	Name	Price	+
124364	Illuminator Plastic White	<b>US\$ 126.38</b>	<a href="#">+Info</a>
124365	Illuminator Gray Enam finish	<b>US\$ 179.52</b>	<a href="#">+Info</a>

**X-ray viewer - TPC**

X-ray viewer - TPC



**Accessories**

Code	Name	Price	+
469760	TPC - FLAT BRACKET FOR UNIT MOUNT (FOR X-RAY VIEWER) - #XV150-F	<b>US\$ 73.60</b>	<a href="#">+Info</a>
469761	TPC - ROUND BRACKET FOR UNIT MOUNT (FOR X-RAY VIEWER) - #XV150-R	<b>US\$ 73.60</b>	<a href="#">+Info</a>

### Other Models

#### Other Models

Code	Name	Price	+
464062	TPC Super Thin LED X-Ray Viewer	US\$ 112.00	<a href="#">+Info</a>
469758	LED X-RAY VIEWER 110V	US\$ 112.00	<a href="#">+Info</a>
469759	LED X-RAY VIEWER (w/Euro Adaptor)	US\$ 112.00	<a href="#">+Info</a>
469762	LED X-RAY VIEWER HOLDER CLIPS 3/PK	US\$ 12.72	<a href="#">+Info</a>

### Parts for X-Ray Viewer TPC

Code	Name	Price	+
469763	TPC - Replacement LED Light Strip (For X-ray Viewer LED Type) - #XV130	US\$ 24.00	<a href="#">+Info</a>

### X-ray viewer - WOLF MEDICAL

#### X-ray viewer - WOLF MEDICAL

Code	Name	Price	+
420843	X-Ray View box Large	US\$ 258.56	<a href="#">+Info</a>

## Catalog Index / Indice

Developer .....	1
Film & Sensors Holders .....	2
Film Hangers .....	8
Film Mounts .....	9
Film Processors .....	12
Films .....	13
Panoramic X-Ray .....	14
Portable X-Rays .....	16
Protectors /Aprons /Ponchos .....	20
X-Ray Accessories .....	25
X-Ray Bite Blocks .....	28
X-Ray Plates Scanners .....	29
X-Ray Sensors .....	32
X-Ray Units .....	35
X-ray viewer .....	36