

Ovens All Types

Heating elements

Heating elements - ECCO

Heating elements - ECCO

Code	Name	Price	+
474345	Heating Element for Sinosteel	US\$ 336.00	+Info

HEATING PLATES FOR OVENS - WHIP MIX

HEATING PLATES FOR OVENS - WHIP MIX

Replacement heating plates assembly for ovens



MEDIUM

Code	Name	Price	+
115058	INFINITY M30 Heating Plate Assembly Side - Set of 2 110vol ts - (Medium) #33915	US\$ 335.59	+Info

LARGE

Code	Name	Price	+
120596	WHIP MIX Heating Plate Assy Rear for L30 - Set Of 2 110v (large) - #33917 (heat...	US\$ 306.54	+Info
120561	WHIP MIX Heating Plate Assy Rear - Set Of 2 110v (large) (12) - #15296	US\$ 406.14	+Info
120547	WHIP MIX Heating Plate Assy Side - Set Of 2 110v (large) (12) - #15297	US\$ 460.38	+Info

Ovens All Types

Muffles for Porcelain Ovens

MUFFLES & HEATING ELEMENTES FOR NEY VULCAN -

All Muffles for Ney Vulcan ovens available.

Features

The difference between the heating element assembly and the muffle assembly is that the heating element assembly is only the heating element plate and nothing else. The Muffle assembly is the entire heating chamber with both heating plates attached.

Make sure you pick the correct voltage for your oven Voltages are grouped

100-110-115-120 volts are all the same

210-220-230 volts are all the same

If you read Bi voltage, means it is the same item for both voltages (set up may be required)



MUFFLES ASSEMBLY

Code	Name	Price	+
103243	NEY-VULCAN - Vulcan 3-130 - Muffle Assy -130 Bivolt 110vol ts 220volts The Muffl...	US\$ 1,169.84	+Info
103245	NEY-VULCAN - Vulcan 3-550 - Muffle Assy -550 110vol ts (20 in x16 in x16 in= 19l...	US\$ 947.54	+Info
103247	NEY-VULCAN - Vulcan 3-550 - Muffle Assy -550 230-220volts (20 in x16 in x16 in= ...	US\$ 1,224.79	+Info
103249	NEY-VULCAN - Vulcan 3-550PD - Muffle Assy 3-550PD- 230-220volts The Muffle assem...	US\$ 1,156.88	+Info
103250	NEY-VULCAN - Vulcan 3-550A - Muffle Assy 3-550A- 230-220volts The Muffle assembl...	US\$ 1,155.46	+Info

MUFFLES - NEY OVENS

MUFFLES - NEY OVENS

We provide quality quartz and non-quartz muffles at a reasonable price.

The first quartz muffle to be offered to the dental industry was for the Ney 660 furnace.

We purchase through the original manufacturers

The ones that are discontinued, we buy them from a distinctive manufacturer.

Our aftermarket Muffle provider, since the Ney 660 was launched, has established a niche among the large furnace manufacturers by providing quality quartz and non-quartz Muffles for most porcelain furnaces on the market today at unbeatable prices far below the OEM and all with outstanding service.



NEY CENTURION MUFFLES

Code	Name	Price	+
103225	NEY-CENTURION - Centurion VPM - Muffle Quartz Spiral - Window 100/120-110vol ts ...	US\$ 650.02	+Info
103228	NEY-CENTURION - Centurion VPM - Muffle Quartz Spiral - Window 230-220volts #9493698	US\$ 650.02	+Info
103222	NEY-CENTURION - Centurion VPM - Muffle Quartz Sleeve - Window 230-220volts #9493960	US\$ 650.02	+Info
103224	NEY-CENTURION - Centurion VPC - Muffle Quartz Spiral - Window 110vol ts / 110Vol...	US\$ 650.02	+Info
103227	NEY-CENTURION - Centurion VPC - Muffle Quartz Spiral - Window 230-220volts #9493698	US\$ 650.02	+Info
103218	NEY CENTURION - Centurion VPC - Muffle Quartz Sleeve - Window 100/120-110vol ts ...	US\$ 650.02	+Info
103234	NEY-CENTURION - Centurion Q-EX - Muffle for Centurion Q-EX 110 volts /120-volts	US\$ 650.02	+Info
103223	NEY-CENTURION - Centurion Q-100 - Muffle Quartz Spiral - Window 110/110Volts #EV...	US\$ 650.02	+Info
103217	NEY CENTURION - Centurion Q-100 - Muffle Quartz Sleeve - Window 110/110Volts #EV...	US\$ 650.02	+Info
103220	NEY-CENTURION - Centurion Q-100 - Muffle Quartz Sleeve - Window 230/220volts #94...	US\$ 650.02	+Info

NEY CERAMCO FINESSE MUFFLE

Code	Name	Price	+
112453	CERAMCO - Muffle for Finesse 110vol ts	US\$ 978.03	+Info

NEY CERAMPRESS MUFFLES

Code	Name	Price	+
103236	NEY-CERAMPRESS - Muffle CCFII- 120-110vol ts	US\$ 668.35	+Info
103238	NEY-CERAMPRESS - Muffle CCFII- 230-220volts #9494091	US\$ 763.55	+Info

NEY MARK MUFFLES

Code	Name	Price	+
103260	NEY - Mark II - Muffle-for Mark2 - 220volt -	US\$ 611.87	+Info
103262	NEY - Mark IV - Muffle-for Mark4 - 220volt -	US\$ 628.45	+Info

NEY STARFIRE

Code	Name	Price	+
112247	Ney Starfire Muffle for 110vol ts	US\$ 611.88	+Info


NEY SUNFIRE MUFFLES

Code	Name	Price	+
112258	Muffle for Sunfire 40 220v	US\$ 705.61	+Info
115834	Muffle for Sunfire 40 110v	US\$ 705.61	+Info

NEY SYSTEM 8 MUFFLE

Code	Name	Price	+
112865	Ney - Muffle for Ney System 8 - 110vol The quartz muffle for Ney System 8- Sunfi...	US\$ 611.90	+Info

MUFFLES FOR ECCO
MUFFLES - ECCO



Code	Name	Price	+
114466	ECCO - Muffles for New Century Furnace - 110vol ts	US\$ 119.00	+Info

MUFFLES FOR OVENS - EDG
MUFFLES FOR OVENS - EDG



We provide quality quartz and non quartz muffles at a reasonable price. The first quartz muffle to be offered to the dental industry was for the Ney 660 furnace. We purchase through the original manufacturers. The ones that are discontinued, we buy them from a distinctive manufacturer. Our aftermarket Muffle provider, since the Ney 660 was launched, has established a niche among the large furnace manufacturers by providing quality quartz and non-quartz Muffles for most porcelain furnaces on the market today at unbeatable prices far below the OEM and all with outstanding service.

Code	Name	Price	+
101446	EDG FV-PLUS SINTER Quartz Muffle for Porcelain Furnace 110-220volts - w/ Instruc...	US\$ 720.95	+Info
101305	EDG FV-PLUS AUTO Quartz Muffle for Porcelain Furnace 110-220volts - w/ Instructi...	US\$ 720.95	+Info

MUFFLES FOR OVENS - JELRUS
MUFFLES FOR OVENS - JELRUS



Description:
We Provide Quality Quartz And Non Quartz Muffles At A Reasonable Price. The first quartz muffle to be offered to the dental industry was for the Ney 660 furnace. We purchase through the original manufacturers. The ones that are discontinued, we buy them from a distinctive manufacturer. Our aftermarket Muffle provider, since the Ney 660 was launched, has established a niche among the large furnace manufacturers by providing quality quartz and non-quartz Muffles

for most porcelain furnaces on the market today at unbeatable prices far below the OEM and all with outstanding service.

Code	Name	Price	+
112269	Muffe-for Vip Vision Executive 300 and Jelrus Wizard & Vista - 220volts	US\$ 738.67	+Info
112394	Muffe-for Vip Vision Executive 300 and Jelrus Wizard & Vista - 110vol ts #23702	US\$ 990.58	+Info
114070	Muffe-for Jelrus Vista 220vol	US\$ 738.67	+Info
101126	JELRUS - Muffle-for Jelrus Vista 110vol - #23702	US\$ 738.67	+Info
112887	Muffe-for VIP X Press - 220volts - Quartz Original Jelrus If your Furnace uses ...	US\$ 1,191.82	+Info
102465	JELRUS -Muffle-for VIP X Press - 110vol ts Quartz Original Jelrus - If your Furn...	US\$ 1,029.02	+Info
111973	Muffe-for Jelrus Megastar 220volts Mufla	US\$ 738.67	+Info
107-JELMUF	JELRUS - Muffle-for Jelrus Megastar 110vol ts #23702 Mufla	US\$ 738.67	+Info
107-149934	JELRUS - Muffle-for JELRUS Ceramat-S - 110vol ts - also available 220volts upon ...	US\$ 738.67	+Info
107-ALPMU	JELRUS - Muffle-for Jelrus Alpha Plus 110/220volt or 110vol t - (when ordered 2...	US\$ 932.40	+Info
102304	Muffe-for Vip Vision Executive 500 - 110vol ts - 110V-#23702 220V-03	US\$ 593.46	+Info
112378	Muffe-for Vip Vision Executive 500 - 220volts -	US\$ 812.52	+Info

MUFFLES FOR OVENS - WHIP MIX

MUFFLES FOR OVENS - WHIP MIX



The first quartz muffle to be offered to the dental industry was for the Ney 660 furnace. We purchase through the original manufacturers

The ones that are discontinued, we buy them from a distinctive manufacturer.

Our aftermarket Muffle provider, since the Ney 660 was launched, has established a niche among the large furnace manufacturers by providing quality quartz and non-quartz Muffles for most porcelain furnaces on the market today at unbeatable prices far below the OEM and all with outstanding service.

FIRELITE

Code	Name	Price	+
120545	WHIP MIX Original #10003 Muffle Solid State Relay - #96506	US\$ 444.63	+Info

INSULATION MUFFLE

Code	Name	Price	+
120526	WHIP MIX Original #10219 Muffle Insulation Side Inner - #96519	US\$ 774.40	+Info
120542	WHIP MIX Original #10214 Muffle Insulation Inner Top - #96514	US\$ 575.69	+Info
120544	WHIP MIX Original #10216 Muffle Insulation Bottom Inner - #96516	US\$ 575.69	+Info
120538	WHIP MIX Original #10217 Muffle Insulation Back Inner - #96517	US\$ 210.77	+Info
120541	WHIP MIX Original #10213 Muffle Insulation Outer Top - #96513	US\$ 160.80	+Info
120543	WHIP MIX Original #10215 Muffle Insulation Bottom Outer - #96515	US\$ 147.65	+Info

PRO MUFFLES

Code	Name	Price	+
120484	WHIP MIX Original #01326 Muffle Pro 200/100 Plus/100 Aft03/06 110v - #96035	US\$ 1,378.38	+Info
110735	WHIP MIX muffle for Pro 100 - 110vol ts #96035	US\$ 1,378.38	+Info
110742	WHIP MIX #01326 Muffle PRO 200/100 PLUS/100 AFT03/06 110vol #96035	US\$ 1,378.38	+Info

Muffles for Porcelain Ovens - AMANN GIRRBACH

Muffles for Porcelain Ovens - AMANN GIRRBACH

Code	Name	Price	+
467318	Ceramill Argovent 2	US\$ 3,004.56	+Info

Muffles for Porcelain Ovens - WHIP MIX

Muffles for Porcelain Ovens - WHIP MIX

Code	Name	Price	+
421022	FIRELITE Muffle Only 220v	US\$ 988.42	+Info

MUFFLES VIP OVENS - IMAGINATION

MUFFLES VIP OVENS - IMAGINATION



Code	Name	Price	+
119571	IMAGINATION Muffle-for IMAGE - 110vol ts - Quartz Original A4159	US\$ 861.03	+Info
119572	IMAGINATION Muffle-for IMAGE - 220volts - Quartz Original - A4159	US\$ 861.01	+Info
119570	IMAGINATION Muffle-for IMAGE XPRESS - 110vol ts - Quartz Original Jelrus If your...	US\$ 892.50	+Info
119569	IMAGINATION Muffle-for IMAGE XPRESS - 220volts - Quartz Original Jelrus If your ...	US\$ 892.50	+Info

Ovens All Types
Parts for ovens

BURN OUT TRAYS - ADS
BURNS OUT TRAYS - ADS

Features

- * Withstands repeated thermal cycles.
- * Will not fracture in furnace.
- * Specially designed channels drain off and hold melted wax and allow hot gasses to escape efficiently.
- * Porous surface provides air circulation and guarantees complete wax combustion.
- * Burn Out Trays can remain in the furnace and protect the furnace muffle from damage.
- * Compatible with all modern muffle materials.
- * Will not block heat transfer.



Code	Name	Price	+
117699	ADS Burn Out Trays 6in X 8in pk (2) - #B304-5	US\$ 57.92	+Info
117698	ADS Burn Out Trays 5in x 6in pk (2) - #B304-4	US\$ 51.23	+Info
117697	ADS Burn Out Trays 5x5 pk (2) - #B304-3	US\$ 46.62	+Info
117696	ADS Burn Out Trays 4in x 5in pk (2) - #B304-2	US\$ 45.70	+Info
117695	ADS Burn Out Trays 4in x 4in pk (2) - #B304-1	US\$ 47.43	+Info

NEY VULCAN REPLACEMENT PARTS
VULCAN - REPLACEMENT PARTS



SWITCHES

Code	Name	Price	+
103027	VULCAN - Switch Pushbutton 1 Led - #9306019	US\$ 5.12	+Info

TRAYS

Code	Name	Price	+
102969	VULCAN - BOTTOM TRAY (floor) - 550 - # 9353057	US\$ 224.21	+Info
102968	VULCAN - BOTTOM TRAY (floor) - 130 - # 9353053	US\$ 161.82	+Info

TUBES

Code	Name	Price	+
102971	VULCAN - Tube- Support- 1750 Shelf - #9353071	US\$ 12.13	+Info

Parts for ovens - ADS

Parts for ovens - ADS

Code	Name	Price	+
118076	Green Guard Pst 2oz	US\$ 15.78	+Info

Parts for ovens - AMANN GIRRBACH

Parts for ovens - AMANN GIRRBACH

Code	Name	Price	+
116144	Heating element Ceramill Therm	US\$ 606.42	+Info
421023	Spindle collet 3mm for Ceramill Mot	US\$ 302.61	+Info
466303	Heating Elements for T3 & AT2 4pk	US\$ 1,572.94	+Info

Parts for ovens - BEGO

Parts for ovens - BEGO

Code	Name	Price	+
418958	Switch Coupler w/ spring Miditherm	US\$ 74.28	+Info
418960	Excenter door switch Miditherm	US\$ 74.28	+Info

Parts for ovens - BUFFALO

Parts for ovens - BUFFALO

Code	Name	Price	+
424560	Hot Plate, UNC Speci	US\$ 14.83	+Info

Parts for ovens - ECCO

Parts for ovens - ECCO

Code	Name	Price	+
474318	Elity 70 Trigger	US\$ 121.00	+Info
474904	Firing Platform PLATE	US\$ 0.02	+Info

Parts for ovens - GENERIC

Parts for ovens - GENERIC

Code	Name	Price	+
470297	DC Control Solid State	US\$ 72.00	+Info
470298	DC Control Solid State Relay	US\$ 72.00	+Info

Parts for ovens - IBEX

Parts for ovens - IBEX

Code	Name	Price	+
469890	Ibex Summit door brick	US\$ 91.80	+Info

469891	lbex Summit memory stick	US\$ 19.20	+Info
469892	lbex heat tray, black anodized.	US\$ 41.97	+Info

Parts for ovens - JELRUS

Parts for ovens - JELRUS

Code	Name	Price	+
470927	Silver Calibration Wire	US\$ 68.25	+Info
473681	VIP-Universal XPress KeyBoard	US\$ 852.60	+Info

Parts for ovens - NABERTHERM

Parts for ovens - NABERTHERM

Code	Name	Price	+
420402	Heating Element for Nabertherm LHT	US\$ 1,326.02	+Info

Parts for ovens - SHENPAZ

Parts for ovens - SHENPAZ

Code	Name	Price	+
464885	Shenpaz Sintra Heating Element	US\$ 650.02	+Info
470223	Sintering table	US\$ 616.00	+Info
474495	Sintra Ceramic Table	US\$ 165.00	+Info
474496	Assembled FOCUS platform	US\$ 295.01	+Info
474848	Sintra Plus Calibration disks-15pk	US\$ 120.00	+Info
475061	Shenpaz Sintra PCB	US\$ 703.80	+Info
476661	Shenpaz Sintra PCB	US\$ 956.80	+Info
476665	SD Card to SINTRA 90	US\$ 210.00	+Info

Parts for ovens - WHIP MIX

Parts for ovens - WHIP MIX

Code	Name	Price	+
120476	#11642 Overlay Pro Press Sp	US\$ 134.00	+Info
120494	#11480 Logic Board Pro Press S	US\$ 1,242.20	+Info
120510	#01561 Pro Press Sp Regulator	US\$ 295.87	+Info
476480	Whipmix pro press 100 thermocouple	US\$ 558.82	+Info

PCB REPLACEMENT PARTS - WHIP MIX

PCB REPLACEMENT PARTS - WHIP MIX

Replacement for Infinity Oven and PDQ Digital Oven



Code	Name	Price	+
120560	WHIP MIX Pcb Replacement Infinity 110v (medium) - #15500-041	US\$ 1,114.37	+Info
118788	WHIP MIX Controller PCB Infinity M-30 110vol #15500-041	US\$ 1,114.37	+Info
120554	WHIP MIX Pcb Replacement- Infinity 110v (large) (12) - #15305	US\$ 902.81	+Info
120556	WHIP MIX Pcb Replacement- Pdq Digital 110v (large) (12) - #15300	US\$ 765.26	+Info
120559	WHIP MIX Pcb Replacement- Pdq Digital 110v (medium) (12) - #15298	US\$ 765.26	+Info
120555	WHIP MIX Pcb Replacement- Infinity 110v (medium) (12) - #15303	US\$ 902.81	+Info
120550	WHIP MIX #33920 Overlay Infinity Large (12) - #15310	US\$ 312.74	+Info
120551	WHIP MIX #33919 Overlay Infinity Medium (12) - #15309	US\$ 257.00	+Info

PORCELAIN FIRING PLATFORM FOR FURNACE - IMAGE

PORCELAIN FIRING PLATFORM FOR FURNACE - IMAGE

Porcelain Firing Platform for Image Press Furnace



Code	Name	Price	+
115919	IMAGE- X-Press Porcelain Firing Press Platform for Image Press- F022 - (Pressin...	US\$ 64.00	+Info
113524	Porcelain Firing Platform F024 -	US\$ 64.00	+Info

PRO 200 REPLACEMENT PARTS - WHIP MIX

PRO 200 REPLACEMENT PARTS - WHIP MIX

These are replacement parts for Pro 200 Serie Furnances



Code	Name	Price	+
120515	WHIP MIX #10931 Solid State Relay 25 Amp - #96225	US\$ 194.23	+Info
120471	WHIP MIX #10343 Firing Tray Regular for PRO - #96015	US\$ 150.39	+Info
120511	WHIP MIX #10360 Cooling Tray - #96003	US\$ 76.59	+Info
120501	WHIP MIX #5560 Lift Limit Sensor - #96094	US\$ 58.41	+Info

120509	WHIP MIX #10369 Power Cord Pump 110v - #96006	US\$ 46.67	+Info
120506	WHIP MIX #11549 250 Volt / 7 Amp Slo-blo Fuse (2pk) - #96092	US\$ 28.72	+Info

PRO REPLACEMENT PARTS - WHIP MIX

PRO REPLACEMENT PARTS - WHIP MIX

These are replacement parts for Pro Serie Furnances



Code	Name	Price	+
120491	WHIP MIX PCB #11481 Logic Board Pro 200 - #96342	US\$ 1,392.24	+Info
120488	WHIP MIX #01009 Vac Valve Assy Pro 200 - #96069	US\$ 569.22	+Info
120492	WHIP MIX #01543 Vacuum Sensor Assembly Pro 200 Series - #96341	US\$ 180.37	+Info
118768	WHIP MIX #11482 KEYPAD ASSEMBLY PRO 200 #96636	US\$ 141.94	+Info
120478	WHIP MIX #11482 Keypad Assembly Pro 200 Series - #96336	US\$ 141.94	+Info
120508	WHIP MIX Flash Drive Pro 200 - #96091	US\$ 132.18	+Info

PRO-PRESS REPLACEMENT PARTS - WHIP MIX

PRO-PRESS REPLACEMENTS PARTS - WHIP MIX

These spare parts will make their dental units more efficient by replacing the broken parts or that do not work for new parties to extend the life of use of your machine



Code	Name	Price	+
120503	WHIP MIX PCB #11481 Logic Mother Board Pro Press 200 - #96343	US\$ 1,392.24	+Info
120507	WHIP MIX #10370 Firing Tray Press w/o Insert - #96013	US\$ 156.75	+Info
120480	WHIP MIX #10371 Insert F/Firing Tray Press for PRO-Press - #96012	US\$ 137.01	+Info
120498	WHIP MIX Flash Drive Pro Press 200 - #96098	US\$ 132.18	+Info

SINTER OVEN REPLACEMENT PARTS - ECCO

SINTER OVEN REPLACEMENT PARTS - ECCO

These spare parts will make their dental units more efficient by replacing the broken parts or that do not work for new parties to extend the life of use of your machine



ECCO

ECCO

Code	Name	Price	+
113539	Firing Platform Stand	US\$ 270.40	+Info
114076	Special Cable #2	US\$ 8.77	+Info
114078	Door Tray for Ecco Zirc	US\$ 48.00	+Info
114952	Platform Pegs for Ecco Zirc 6pk	US\$ 15.36	+Info
114953	Sockets for heating elements	US\$ 19.20	+Info
114954	5 Amp Fuse for Ecco-Zirc	US\$ 4.80	+Info
114955	32 Amp Fuse for Ecco-Zirc	US\$ 4.80	+Info
114956	Platform & Fans Circuit board	US\$ 80.00	+Info
114974	Relay for Sinter Oven	US\$ 19.20	+Info
114975	Aluminium Connection Cable #1	US\$ 8.77	+Info
115669	Zirconia Sinter Oven Motor	US\$ 320.00	+Info
115912	Main Mother Circuit board 220v	US\$ 368.00	+Info
116678	Muffle Chamber for Ecco-Sint	US\$ 3,067.93	+Info

SINOSTEEL


SINOSTEEL

Code	Name	Price	+
113796	Heating Element for Ecco Zirc	US\$ 240.00	+Info

REPLACEMENT PARTS FOR 500S FURNANCE - DLS

REPLACEMENT PARTS FOR 500S FURNANCE - DLS

These spare parts will make their dental units more efficient by replacing the broken parts or that do not work for new parties to extend the life of use of your machine



Code	Name	Price	+
110658	DLS - Control Board for Model 500S Furnace - Burnout Furnace #E520 Replacement P...	US\$ 336.00	+Info
110656	DLS - Heating Elements set for Model 500S Furnace #E506 Replacement Parts for 50...	US\$ 247.50	+Info
110664	DLS - Door Assembly (with packing) for Model 500S Furnace - Burnout Furnace #E5...	US\$ 189.42	+Info
110659	DLS - Triac Assembly for Model 500S Furnace - Burnout Furnace #E522 Replacement ...	US\$ 102.90	+Info
110660	DLS - Classic Triac Assembly for Model 500S Furnace - Burnout Furnace #E524 Repl...	US\$ 82.95	+Info
110661	DLS - Triac Only for Model 500S Furnace - Burnout Furnace #E526 Replacement Part...	US\$ 68.25	+Info

REPLACEMENTS PARTS FOR OVENS - ECCO

REPLACEMENTS PARTS FOR OVENS - ECCO

These spare parts will make their dental units more efficient by replacing the broken parts or that do not work for new parties to extend the life of use of your machine



Code	Name	Price	+
118765	ECCO Transformer 220 volsts for 3F ECCO sinterization oven	US\$ 480.00	+Info
118766	ECCO - Main Trigger Board for 3F ECCO-Zirc Sinterization Furnace 220volts - each-	US\$ 576.00	+Info
102502	ECCO - Digital Temp Controler for Digital CFS furnace fort 110vol ts 220 v (Fle...	US\$ 70.40	+Info
114560	ECCO - Firing Platform New Century	US\$ 99.09	+Info
122323	ECCO - Zirconia Firing Platform SWITCH For ECCO Sinterilization Furnace -	US\$ 32.00	+Info
117599	ECCO - Sinter Furnace Fan each (Parts and accesories for ECCO - Sinter Furnace)	US\$ 24.00	+Info

REPLACEMENTS PARTS FOR OVENS - JELRUS

REPLACEMENTS PARTS FOR OVENS - JELRUS

These spare parts will make their dental units more efficient by replacing the broken parts or that do not work for new parties to extend the life of use of your machine



Code	Name	Price	+
107-ALPTER	JELRUS - Thermocouple for Alpha Plus 220 volts - #Y061	US\$ 493.50	+Info
110731	JELRUS - Vip X Press Press Ring insert Kit - #A4153	US\$ 362.26	+Info
110730	JELRUS - Vip X Press Press Firing Block Kit- #A4152	US\$ 266.82	+Info
110629	JELRUS - Firing Platform Jelrus Vista - NOW DISCONTINUED - A NEW CENTURY Furnace...	US\$ 104.00	+Info
107-JELORI	JELRUS - Furnance O-Rings Gasket Door/Muffle Top 46 inches Ea. #36039 Anillo Sel...	US\$ 27.22	+Info
124120	JELRUS - VIP-300 VISION KeyPad Panel Membrane #24525	US\$ 115.20	+Info
123517	Jelrus Image Press PCBoard 110 #A4178	US\$ 1,217.97	+Info
123518	Jelrus Image Press Doors Oinrings #A4182	US\$ 72.45	+Info

REPLACEMENTS PARTS FOR OVENS - WHIP MIX

REPLACEMENTS PARTS FOR OVENS - WHIP MIX

These spare parts will make their dental units more efficient by replacing the broken parts or that do not work for new parties to extend the life of use of your machine



DOOR REPLACEMENT PARTS

Code	Name	Price	+
120585	WHIP MIX Door Hinge Kit - Set Of 2 - #33998	US\$ 42.72	+Info
120573	WHIP MIX Door Hook Spring Assy Kit - #33997	US\$ 48.17	+Info
120563	WHIP MIX #33922 Door Switch (12) - #15294	US\$ 63.20	+Info
120570	WHIP MIX Door Handle Kit (medium) - #15279	US\$ 128.51	+Info
120592	WHIP MIX Door Insulation (medium) - #15722	US\$ 265.80	+Info
120593	WHIP MIX Door Insulation (large) - #15712	US\$ 298.49	+Info
120575	WHIP MIX Door Assembly (medium) - #15284	US\$ 537.15	+Info
120567	WHIP MIX Door Assembly (large) - #15286	US\$ 562.93	+Info
120536	WHIP MIX #10212 Door Insulation - #96512	US\$ 982.45	+Info

SWITCH REPLACEMENT PARTS

Code	Name	Price	+
120481	WHIP MIX #11246 Limit Switch - #96026	US\$ 58.38	+Info
120594	WHIP MIX On/off Power Switch 110v / 220vol - #15684	US\$ 41.58	+Info

CERAMIC FRONT SECTION

Code	Name	Price	+
120564	WHIP MIX Ceramic Front Section - Set Of 2 (large) - #15293	US\$ 143.77	+Info
120565	WHIP MIX Ceramic Front Section - Set Of 2 (medium) - #15292	US\$ 107.27	+Info

FLOOR PLATES

Code	Name	Price	+
120548	WHIP MIX Floor Plate w/filler Strip Insulation (large) - #33981	US\$ 156.08	+Info

HEATER REPLACEMENT PARTS

Code	Name	Price	+
120583	WHIP MIX #33929 Heater Wire Assembly (12) - #15312	US\$ 93.83	+Info
120589	WHIP MIX Heater-jumper - #33130	US\$ 13.40	+Info

INSULATING REPLACEMENT PARTS

Code	Name	Price	+
120576	WHIP MIX Ceramic Insulating Bushing - Set Of 4 - #33958	US\$ 19.99	+Info

Other Models

Other Models

Code	Name	Price	+
120475	Flash Drive Pro Press Sp	US\$ 125.74	+Info
120495	#10239 Lift Table	US\$ 382.35	+Info

OTHER REPLACEMENTS PARTS

Code	Name	Price	+
120539	WHIP MIX PCB #33907 Logic Board - #96504	US\$ 439.73	+Info
120557	WHIP MIX Circuit Breaker - 110v 16 Amp - #117046	US\$ 101.64	+Info
120587	WHIP MIX Tray For Heating Chamber - #33256	US\$ 66.21	+Info
120566	WHIP MIX Calibration Table Kit 1500f (816c) pk 5 - #15291	US\$ 27.45	+Info
120591	WHIP MIX Vent Tube - #15729	US\$ 27.16	+Info

Ovens All Types

Thermocouples



REPLACEMENT THERMOCOUPLE S FOR OVENS - WHIP MIX



Code	Name	Price	+
116082	WHIP MIX Thermocouple for Burnout Furnace Infinity M-30 L-30 Thermocouple insulated	US\$ 169.50	+Info
120590	WHIP MIX Thermocouple Assembly Kit - #18913	US\$ 169.50	+Info
120562	WHIP MIX #33963 Thermocouple Assy With Connector (12) - #15295	US\$ 238.81	+Info
120552	WHIP MIX #33918 Overlay Pdq Large (12) - #15308	US\$ 307.64	+Info
118767	WHIP MIX Thermocouple for Pro Press 200 insulated #96312	US\$ 914.90	+Info
120528	WHIP MIX #10001 Heating Elements pk 6 - #96510	US\$ 3,094.18	+Info

Thermocouple for BEGO MIDTHERM



Code	Name	Price	+
10-14087	BEGO Thermocouple Miditherm 100/200 MP 2pk #14087	US\$ 317.30	+Info

THERMOCOUPLE FOR VULCAN - NEY

THERMOCOUPLE FOR VULCAN - NEY



Code	Name	Price	+
103206	NEY VULCAN - Vulcan -1750 - Thermocouple Assembly - Series II & 1750 A-1750 3-17...	US\$ 247.00	+Info

Thermocouples - DLS

Thermocouples - DLS

Code	Name	Price	+
110657	Thermocouple for Model 500S	US\$ 91.29	+Info

Thermocouples - ECCO

Thermocouples - ECCO

Code	Name	Price	+
114806	Thermocouple for New Century	US\$ 324.14	+Info

Thermocouples - EDG

Thermocouples - EDG

Code	Name	Price	+
30699	Alumine 50 Thermocouple	US\$ 90.84	+Info

Thermocouples - JELRUS

Thermocouples - JELRUS

Code	Name	Price	+
100825	Thermocouple Jelrus Infinity	US\$ 166.33	+Info
113802	Thermocouple Alpha Plus 110v	US\$ 308.66	+Info
114957	Circuit Board Jelrus Inf M-30	US\$ 760.36	+Info
122998	Thermocouple Jelrus VIP EXPRES	US\$ 573.30	+Info
473699	Thermocouple Radiance Multistage	US\$ 166.33	+Info
473816	Thermoc Jelrus Wizard/Vision/Exec	US\$ 576.24	+Info
474512	Thermocouple Shenpaz 4010 porc	US\$ 280.01	+Info

Thermocouples - SINOSTEEL

Thermocouples - SINOSTEEL

Code	Name	Price	+
113797	Thermocouple for Ecco Zirc	US\$ 720.00	+Info

Catalog Index / Indice

Heating elements	1
Muffles for Porcelain Ovens	1
Parts for ovens	6
Thermocouples	15