

Abutments HEX External Hexagon

Abutments External Hex



ABUTMENTS EXTERNAL HEX - XPS

ABUTMENT EXTERNAL HEX (Branemark type)

* **Platform 3.5mm** is compatible with following Brands 3i Nanotite & Osseotite BRANEMARK B&W ,IMPLANT FORT & TREE OSS

* **Platform 4.1mm** is compatible with following Brands : 3i Nanotite & Osseotite , BIOHORIZONS , BRANEMARK , B&W IMPLANT FORT, TREE OSS , LIFERCORE , Steri-Oss HL , Center Pulse Taper-Lock

* **Platform 5.0mm** is Compatible with following Brands and Diameters :
 3i Nanotite & Osseotite 5.0 Diam / 5.0 Platform
 BRANEMARK 3.75 & 4.0 Diam / 4.1 Platform
 B&W 5.0 Diam / 5.0 Platform
 TREE OSS 5.0 Diam / 5.1 Platform



NEEDED INSTRUMENTS

Same instruments than used for NR
 Wrench (Square engagement)
 Fixation Screw 1.2mm (0.48in)
 Screwdriver 0.048in 1.2mm for RP and WP and
 Screwdriver 0.035in 0.9mm for NP

Real shape might differ from the diagram. To check real shapes please click on each individual item.

Branemark and Nobel are a Trademark of Biocare.

TITANIUM STRAIGHT Abutment

3.5mm

Code	Name	Price	+
110751	XPS External Hex Titanium Straight Abutment 3.5 Platform / 2mm Collar Height use...	US\$ 18.00	+Info

5.0mm

XGATE Dental

XGATE Dental

Code	Name	Price	+
------	------	-------	---

110752	Ext Hex Tit Str Abut 5.0/1mm	US\$ 18.00	+Info
--------	------------------------------	------------	-----------------------

XPS
XPS

Code	Name	Price	+
110753	Ext Hex Tit Str Abut 5.0/2mm	US\$ 18.00	+Info

TITANIUM ANGLED Abutment

15

Code	Name	Price	+
110786	XPS External Hex Titanium Angled Abutment 15 / 4.1mm Platform / 1mm Collar Heig...	US\$ 18.00	+Info

MULTI UNIT STRAIGHT Abutment 4.1mm

Code	Name	Price	+
110811	XPS External Hex Multi-Unit-Platform 4.1 Platform / 1mm-Exp. Profile To 4.8mm Co...	US\$ 18.00	+Info
110812	XPS External Hex Multi-Unit-Platform 4.1 Platform / 2mm-Exp. Profile To 4.8mm Co...	US\$ 18.00	+Info
110813	XPS External Hex Multi-Unit-Platform 4.1 Platform / 3mm-Exp. Profile To 4.8mm Co...	US\$ 18.00	+Info
110814	XPS External Hex Multi-Unit-Platform 4.1 Platform / 4mm-Exp. Profile To 4.8mm Co...	US\$ 18.00	+Info

Other Models
Other Models

Code	Name	Price	+
110750	Ext Hex Tit Str Abut 3.5 /1mm	US\$ 18.00	+Info
110754	Ext Hex Tit Est Abut 5.0/1mm	US\$ 18.00	+Info

Abutments HEX External Hexagon

Analogs External Hex



ANALOGS EXTERNAL HEX - BRANEMARK

LABORATORY ANALOG - EXTERNAL (BRANEMARK)

- * BRONZE Laboratory Analog
- * STEEL Laboratory Analog
- * TITANIUM Laboratory Analog

NOTE

Actual shapes may differ from the schematic drawings shown above.
To check real pictures please click ...Read More on desired item



TITANIUM Laboratory Analog

Code	Name	Price	+
110757	XPS External Hex Titanium Laboratory Analog 3.3mm Platform / (Hex 2.7mm compati...	US\$ 10.00	+Info
110792	XPS External Hex Titanium Laboratory Analog 4.1mm Platform (Hex 2.7mm compatibl...	US\$ 10.00	+Info

110758	XPS External Hex Titanium Laboratory Analog 5.0mm Platform /	US\$ 10.00	+Info
--------	--	------------	-----------------------

Abutments HEX External Hexagon

Fixation Screws External Hex



Fixation Screws External Hex - XPS

Fixation Screws External Hex - XPS

Code	Name	Price	+
110816	Ext Hex Tit Hx Screw RP & WP	US\$ 7.80	+Info

FIXATIONS SCREWS EXTERNAL HEX - XPS

FIXATIONS SCREWS EXTERNAL HEX - XPS

We have 2 versions :

NP (Narrow Platform)

RP (Regular Platform) and WP (Wide Platform)

NP:

Thread Diameter : 1.80x 0.35

Material : Titanium 5

Head Diameter : 2.45mm

Head Length : 2mm

Total Length : 9mm

Wrench type : HEX 0.50 in 1.25mm

RP & WP:

Thread Diameter : 2.00x 0.40

Material : Titanium 4

Head Diameter : 2.45mm

Head Length : 2mm

Total Length : 10.3mm

Wrench type : HEX 0.50 in 1.25mm

All Screws use the same Wrench and are 100% compatible with Nobel Replace screws.



SCREW HEX 0.48

Code	Name	Price	+
110815	XPS External Hex Bronze Laboratory Fixation Screw only, hexagon fit 0.048 in 1pk	US\$ 7.80	+Info

SCREW SQUARE 0.50in

Code	Name	Price	+
110817	XPS External Hex Bronze Laboratory Fixation Screw only, Square fit 0.50 in	US\$ 7.80	+Info
110818	XPS External Hex Titanium Fixation Screw only, Square fit 0.50 in	US\$ 6.00	+Info

Abutments HEX External Hexagon

Healing Caps External Hex



HEALING ABUTMENTS EXTERNAL HEX - XPS

Healing Caps - External (Branemark)

Platform 3.5mm is compatible with following Brands and Diameters :

- 3i Nanotite & Osseotite 3.25 Diam / 3.4 Platform
- BRANEMARK 3.25 Diam / 3.5 Platform
- B&W 3.3 Diam / 3.4 Platform
- IMPLANT FORT 3.3 Diam / 3.4 Platform
- TREE OSS 3.3 Diam / 3.5 Platform



Platform 4.1mm is compatible with following Brands and Diameters :

- 3i Nanotite & Osseotite 4.0 Diam / 4.1 Platform
- BIOHORIZONS 4.6 Diam / 4.5 Platform
- BRANEMARK 3.75 & 4.0 Diam / 4.1 Platform
- B&W 3.75 & 4.0mm Diam / 4.1 Platform
- IMPLANT FORT 3.75 & 4.0 Diam / 4.1 Platform
- TREE OSS 3.75 & 4.0 Diam / 4.1 Platform
- LIFERCORE 3.75, 4.0 Diam / 4.1 Platform
- Steri-Oss HL - 3.8, 4.5 Diam
- Center Pulse Taper-Lock - 3.75 Diam / 4.1 Platform

Platform 5.0mm is Compatible with following Brands and Diameters :

- 3i Nanotite & Osseotite 5.0 Diam / 5.0 Platform
- BRANEMARK 3.75 & 4.0 Diam / 4.1 Platform
- B&W 5.0 Diam / 5.0 Platform
- TREE OSS 5.0 Diam / 5.1 Platform

Expanding Profile Healing Caps available :

- * To 4.5 mm
- * To 5.0 mm
- * To 6.0 mm

NEEDED INSTRUMENTS

- Same instruments than used for NR
- Wrench (Square engagement)
- Fixation Screw 1.2mm (0.48in)
- Screwdriver 0.048in 1.2mm for RP and WP and
- Screwdriver 0.035in 0.9mm for NP

3.5mm PLATFORM

Code	Name	Price	+
110764	XPS External Hex Healing Abutment 3.3mm Platform / 2mm Collar Height - 2.4mm eng...	US\$ 7.90	+Info
110765	XPS External Hex Healing Abutment 3.3mm Platform / 4mm Collar Height - 2.4mm eng...	US\$ 7.90	+Info

4.1mm PLATFORM

Code	Name	Price	+
110799	XPS External Hex Healing Abutment 4.1 Platform / 2mm Collar Height - 2.7mm engag...	US\$ 7.90	+Info
110800	XPS External Hex Healing Abutment 4.1Platform / 3mm Collar Height - 2.7mm engagi...	US\$ 7.90	+Info

110801	XPS External Hex Healing Abutment 4.1Platform / 4mm Collar Height - 2.7mm engagi...	US\$ 7.90	+Info
--------	---	-----------	-----------------------

Abutments HEX External Hexagon
Transfers External Hex



TRANSFERS EXTERNAL HEX BRANEMARK

TRANSFER COPING POST EXTERNAL (BRANEMARK)

Platform 3.5mm is compatible with following Brands and Diameters:

- 3i Nanotite & Osseotite 3.25 Diam / 3.4 Platform
- BRANEMARK 3.25 Diam / 3.5 Platform
- B&W 3.3 Diam / 3.4 Platform
- IMPLANT FORT 3.3 Diam / 3.4 Platform
- TREE OSS 3.3 Diam / 3.5 Platform



Platform 4.1mm is compatible with following Brands and Diameters:

- 3i Nanotite & Osseotite 4.0 Diam / 4.1 Platform
- BIOHORIZONS 4.6 Diam / 4.5 Platform
- BRANEMARK 3.75 & 4.0 Diam / 4.1 Platform
- B&W 3.75 & 4.0mm Diam / 4.1 Platform
- IMPLANT FORT 3.75 & 4.0 Diam / 4.1 Platform
- TREE OSS 3.75 & 4.0 Diam / 4.1 Platform
- LIFERCORE 3.75, 4.0 Diam / 4.1 Platform
- Steri-Oss HL - 3.8, 4.5 Diam
- Center Pulse Taper-Lock - 3.75 Diam / 4.1 Platform

Platform 5.0mm is Compatible with following Brands and Diameters:

- 3i Nanotite & Osseotite 5.0 Diam / 5.0 Platform
- BRANEMARK 3.75 & 4.0 Diam / 4.1 Platform
- B&W 5.0 Diam / 5.0 Platform
- TREE OSS 5.0 Diam / 5.1 Platform
- Platform 5.0mm is Compatible with 3i.

NEEDED INSTRUMENTS

- Fixation Screw 1.2mm (0.48in)
- Screwdriver 0.048in 1.2mm for RP and WP and
- Screwdriver 0.035in 0.9mm for NP
- Branemark is a Trademark of Biocare.

NOTE

*Actual shapes may differ from the schematic drawings shown above.
To check real pictures please click ...Read More on desired item*

[Transfers Closed Tray](#)

XGATE Dental
XGATE Dental

Code	Name	Price	+
110766	Ext Hex Trans Post 3.5mm	US\$ 21.42	+Info

XPS
XPS

Code	Name	Price	+
110767	Ext Hex Trans Post 5.0mm	US\$ 12.80	+Info
110821	Ext Hex Transfer 4.1 Close Tra	US\$ 12.80	+Info

Transfers Open Tray

Code	Name	Price	+
110819	XPS External Hex Impression Transfer Open Tray 4.1 Platform	US\$ 12.80	+Info

Abutments HEX External Hexagon

Ucla External Hex



UCLA 360 EXTERNAL HEX (BRANEMARK)

Material: Acrylic Machined / Machined POM

Innovation: UCLA 360 may be angled manually at any angle between 0 and 20 and after casting, the torque can be performed even after the screw angled.



Use: Used in casting prosthesis directly on the hexagon of the implant securing the prosthesis. Can be used in cemented prostheses or provisional unit for elements (hexagonal) or multiple (rotational).

Castable Plastic Engaging

Code	Name	Price	+
114242	XPS External Hex UCLA-360 Castable No-Rotational w/Crco base 3.4mm Platform / 1m...	US\$ 10.80	+Info
114243	XPS External Hex UCLA-360 Castable No-Rotational w/Crco base 4.1mm Platform / 1m...	US\$ 10.80	+Info
114244	XPS External Hex UCLA-360 Castable No-Rotational w/Crco base 5.0mm Platform / 1m...	US\$ 10.80	+Info

Cr Cobalt Base Engaging

Code	Name	Price	+
114239	XPS External Hex UCLA-360 Castable No-Rotational 3.4mm Platform / 1mm Collar Hei...	US\$ 10.80	+Info
114240	XPS External Hex UCLA-360 Castable No-Rotational 4.1mm Platform / 1mm Collar He...	US\$ 10.80	+Info
114241	XPS External Hex UCLA-360 Castable No-Rotational 5.0mm Platform / 1mm Collar Hei...	US\$ 10.80	+Info

UCLA EXTERNAL HEX (BRANEMARK)

Platform 3.5mm is compatible with following Brands and Diameters:

3i Nanotite & Osseotite 3.25 Diam / 3.4 Platform

BRANEMARK 3.25 Diam / 3.5 Platform

B&W 3.3 Diam / 3.4 Platform

IMPLANT FORT 3.3 Diam / 3.4 Platform

TREE OSS 3.3 Diam / 3.5 Platform



Platform 4.1mm is compatible with following Brands and Diameters:

3i Nanotite & Osseotite 4.0 Diam / 4.1 Platform

BIOHORIZONS 4.6 Diam / 4.5 Platform

BRANEMARK 3.75 & 4.0 Diam / 4.1 Platform

B&W 3.75 & 4.0mm Diam / 4.1 Platform

IMPLANT FORT 3.75 & 4.0 Diam / 4.1 Platform

TREE OSS 3.75 & 4.0 Diam / 4.1 Platform

LIFERCORE 3.75, 4.0 Diam / 4.1 Platform

Steri-Oss HL - 3.8, 4.5 Diam

Center Pulse Taper-Lock - 3.75 Diam / 4.1 Platform

Platform 5.0mm is Compatible with following Brands and Diameters:

3i Nanotite & Osseotite 5.0 Diam / 5.0 Platform

BRANEMARK 3.75 & 4.0 Diam / 4.1 Platform

B&W 5.0 Diam / 5.0 Platform

TREE OSS 5.0 Diam / 5.1 Platform

NEEDED INSTRUMENTS

Same instruments than used for NR

Wrench (Square engagement)

Fixation Screw 1.2mm (0.48in)

Screwdriver 0.048in 1.2mm for RP and WP and

Screwdriver 0.035in 0.9mm for NP

NOTE

Actual shapes may differ from the schematic drawings shown above.

To check real pictures please click ...Read More on desired item

3.5mm PLATFORM

Code	Name	Price	+
110768	XPS External Hex UCLA-Castable Plastic No-Rotational 3.4 Platform / 1mm Collar H...	US\$ 5.80	+Info

4.1mm PLATFORM ENGAGING

Code	Name	Price	+
110824	XPS External Hex UCLA-Castable- Plastic No-Rotational 4.1 Platform / 1mm Collar	US\$ 12.80	+Info
110825	XPS External Hex UCLA-Castable- Plastic No-Rotational 4.1 Platform / 2mm Collar	US\$ 12.80	+Info

4.1mm PLATFORM NO ENGAGING

Code	Name	Price	+
110823	XPS External Hex UCLA-Castable- No-Rotational Platform 4.1 Platform / 2mm Collar	US\$ 12.80	+Info
110822	XPS External Hex UCLA-Castable- yes-Rotational Plastic Platform 4.1 Platform / ...	US\$ 12.80	+Info

5.0mm PLATFORM

Platform 5.0mm is Compatible with 3I.Branemark is a Trademark of Biocare.

Code	Name	Price	+
110769	XPS External Hex UCLA-Castable-Plastic No-Rotational 5.0 Platform / 1mm Collar ...	US\$ 11.00	+Info

UCLA WITH CHROME COBALT BASE



Code	Name	Price	+
110770	XPS External Hex UCLA-B yes-rotational w/Chrom Cobalt Base 3.4mm Platform / 1mm...	US\$ 16.00	+Info
116106	XPS External Hex UCLA-B No-Rotational w/Crco base - 4.1mm Platform / 1mm Collar ...	US\$ 16.00	+Info
110771	XPS External Hex UCLA-B No-Rotational w/Chrom Cobalt Base 5.0mm Platform / 1mm C...	US\$ 16.00	+Info

Catalog Index / Indice

Abutments External Hex	1
Analogs External Hex	2
Fixation Screws External Hex	3
Healing Caps External Hex	4
Transfers External Hex	5
Ucla External Hex	6