

Abutments Internal Hex

Abutments Internal Hex



Abutments Internal Hex - XGATE DENTAL

Abutment Internal Hex - XGATE DENTAL



15 Angled

15 Angled

Code	Name	Price	+
463687	XGATE - IntHex 15 Ang d abutme	US\$ 30.03	+Info
463688	XGATE - IntHex 15 Ang d abutme	US\$ 30.03	+Info
463689	XGATE - IntHex 15 Ang d anatomi	US\$ 30.03	+Info
463690	XGATE - IntHex 15 Ang d anatomi	US\$ 30.03	+Info
463691	XGATE - IntHex 15 Ang d anatomi	US\$ 30.03	+Info
463692	XGATE - IntHex 15 Ang d anatomi	US\$ 30.03	+Info

25 Angled

25 Angled

Code	Name	Price	+
463693	XGATE - IntHex 25 Ang d abutme	US\$ 30.03	+Info
463694	XGATE - IntHex 25 Ang d abutme	US\$ 30.03	+Info
463695	XGATE - IntHex 25 Ang d anatomi	US\$ 30.03	+Info
463696	XGATE - IntHex 25 Ang d anatomi	US\$ 30.03	+Info
463697	XGATE - IntHex 25 Ang d anatomi	US\$ 30.03	+Info
463698	XGATE - IntHex 25 Ang d anatomi	US\$ 30.03	+Info

Anatomic Abutment

Anatomic Abutment

Code	Name	Price	+
463701	IntHex Straight anatomic 1mm	US\$ 18.59	+Info

463702	IntHex Straight anatomic 2mm	US\$ 18.59	+Info
463703	IntHex Straight anatomic 3mm	US\$ 18.59	+Info
463704	IntHex Straight anatomic 4mm	US\$ 18.59	+Info

Regular Straight

Regular Straight

Code	Name	Price	+
463699	IntHex Straight abutment - 11mm	US\$ 14.30	+Info
463700	IntHex Straight abutment 9mm	US\$ 14.30	+Info
463705	IntHex Straight abutment 11mm	US\$ 18.59	+Info
463706	IntHex Straight abutment 9mm	US\$ 18.59	+Info
463707	XGATE - IntHex Straight abutment w	US\$ 18.59	+Info
463708	XGATE - IntHex Straight abutment w	US\$ 18.59	+Info
463709	XGATE - IntHex Straight abutment w	US\$ 18.59	+Info
463710	XGATE - IntHex Straight abutment w	US\$ 18.59	+Info
463715	XGATE - IntHex Straight abutment w	US\$ 18.59	+Info
463716	XGATE - IntHex Straight abutment w	US\$ 18.59	+Info
463717	XGATE - IntHex Straight abutment w	US\$ 18.59	+Info
463718	XGATE - IntHex Straight abutment w	US\$ 18.59	+Info

ANGLED ESTHETIC ABUTMENTS INTHEX - XPS

ANGLED ABUTMENT INTERNAL HEX (XPS)



Compatible with following Brands:

- * MIS seven
- * ZIMMER

Engaging Hexagon 2.45mm (same as Mis & Zimmer)

Platforms:

- * 3.3 expanding to 4.0
- * 4.0 expanding to 4.7
- * 5.0 No expanding
- * Angled : 15 grades

NEEDED INSTRUMENTS

- * Same instruments than used for compatible systems
- * Wrench (Square engagement)
- * Fixation Screw 1.2mm (0.48in)
- * Screwdriver 0.048in 1.2mm for RP and WP and
- * Screwdriver 0.035in 0.9mm for NP

Real shape might differ from the diagram. To check real shapes please click on each individual item.

Code	Name	Price	+
117351	XPS INTHEX Angled Esthetic Abutment 15 for RP 4.0 (Engaging 2.45 - Platform 4...	US\$ 18.00	+Info
117350	XPS INTHEX Angled Esthetic Abutment 15 for NP & RP - (Engaging 2.45 - Platform...	US\$ 18.00	+Info
117352	XPS INTHEX Angled Esthetic Abutment 15 for WP 5.0 - (2.45 mm Engaging - 5.0 Pl...	US\$ 18.00	+Info

STRAIGHT ABUTMENTS INTERNAL HEX INTHEX - XPS

ABUTMENT INTERNAL HEX (XPS)



Compatible with following Brands:

- * MIS seve
- * ZIMMER

Engaging Hexagon 2.45mm (same as Mis & Zimmer)

Platforms:

- 3.3 expanding to 4.0
- 4.0 expanding to 4.7
- 5.0 No expanding

NEEDED INSTRUMENTS

- Same instruments than used for compatible systems
- Wrench (Square engagement)
- Fixation Screw 1.2mm (0.48in)
- Screwdriver 0.048in 1.2mm for RP and WP and
- Screwdriver 0.035in 0.9mm for NP

Real shape might differ from the diagram. To check real shapes please click on each individual item.

Code	Name	Price	+
117339	XPS INTHEX Straight Esthetic Abutment for RP 4.0 (Engaging 2.45 - Platform 4.0 ...	US\$ 15.20	+Info
117340	XPS INTHEX Straight Esthetic Abutment for WP 5.0 (Engaging 2.45 - 5.0 Platform)...	US\$ 15.20	+Info
117341	XPS INTHEX Straight Esthetic Abutment for NP & RP 3.5 Diameter (Engaging 2.45 -...	US\$ 15.20	+Info

Abutments Internal Hex

Analogs Internal Hex



Analogs Internal Hex - XGATE DENTAL

Analogs Internal Hex - XGATE DENTAL



Conical

Conical

Code	Name	Price	+
463747	XGATE - Implant Analog Mini CONICAL	US\$ 18.59	+Info
463750	Imp Anal Reg Conical	US\$ 18.59	+Info

IntHex

IntHex

Code	Name	Price	+
463659	XGATE - Implant analog for Multi U	US\$ 20.02	+Info
463746	XGATE - IntHex Implant analog - An	US\$ 10.01	+Info
463748	XGATE - IntHex Implant analog plus	US\$ 14.28	+Info
463749	XGATE - IntHex Implant analog plus	US\$ 13.99	+Info
463751	XGATE - IntHex Implant analog- wid	US\$ 10.01	+Info
463794	Digitalanalog MUA, D-type	US\$ 38.70	+Info
463795	Digitalanalog MUA, V-type	US\$ 38.70	+Info
463796	Digitalimplant analog	US\$ 38.70	+Info

ANALOGS INTERNAL HEX - XPS

ANALOGS INTERNAL HEX (XPS)

Compatible with following Brands:

- * MIS seven
- * ZIMMER

Engaging Hexagon 2.45mm (same as Mis & Zimmer)

Platforms: 3.3 , 4.0 & 5.0 No expanding

NEEDED INSTRUMENTS

- * Same instruments than used for compatible systems
- * Wrench (Square engagement)
- * Fixation Screw 1.2mm (0.48in)
- * Screwdriver 0.048in 1.2mm for RP and WP and
- * Screwdriver 0.035in 0.9mm for NP

Real shape might differ from the diagram. To check real shapes please click on each individual item.



Code	Name	Price	+
117342	XPS INTHEX Titanium Implant Analog for NP 3.5mm diameter - compatible with MIS &...	US\$ 10.80	+Info
117343	XPS INTHEX Titanium Implant Analog for RP 4.0 - (Engaging 2.45 - Platform 4.0)...	US\$ 10.80	+Info
117344	XPS INTHEX Titanium Implant Analog for WP 5.0 - compatible with MIS & Zimmer	US\$ 10.80	+Info

Abutments Internal Hex

Cad/Cam Abutments

Cad/Cam Abutments - XGATE DENTAL

Cad/Cam Abutments - XGATE DENTAL

Code	Name	Price	+
463793	Scan body implant level	US\$ 42.84	+Info

Abutments Internal Hex
Healing Caps Internal Hex



HEALING CAPS INTERNAL HEX INTHEX - XPS

BALL ABUTMENTS INTERNAL HEX INTHEX - XPS

Compatible with following Brands:

- * MIS seven
- * ZIMMER

Engaging Hexagon 2.45mm (same as Mis & Zimmer)

Platforms:

- 3.3 expanding to 4.0
- 4.0 expanding to 4.7
- 5.0 No expanding



NEEDED INSTRUMENTS

- * Same instruments than used for compatible systems
- * Wrench (Square engagement)
- * Fixation Screw 1.2mm (0.48in)
- * Screwdriver 0.048in 1.2mm for RP and WP and
- * Screwdriver 0.035in 0.9mm for NP

Real shape might differ from the diagram. To check real shapes please click on each individual item.

HEALING ABUTMENT - XPS

Code	Name	Price	+
117334	XPS INTHEX Healing Abutment for NP & RP 3.3 Diameter (Engaging 2.45 - Platform ...	US\$ 7.90	+Info
117335	XPS INTHEX Healing Abutment for RP 4.0 (Engaging 2.45 - Platform 4.0 / Expandin...	US\$ 7.90	+Info
117336	XPS INTHEX Healing Abutment for WP 5.0 Diameter (2.45 mm Engaging - 5.0 Platfor...	US\$ 7.90	+Info

Abutments Internal Hex

Screws



Fixation Screws Internal Hex - XGATE DENTAL

Screw - XGATE DENTAL



Code	Name	Price	+
------	------	-------	---

463776	XGATE - IntHex Screw for abutment	US\$ 3.43	+Info
470508	XGATE - Screw for MU	US\$ 7.14	+Info

FIXATION SCREWS INTERNAL HEX - XPS

We have 2 versions:

NP (Narrow Platform)

RP (Regular Platform) and WP (Wide Platform)

NP:

Thread Diameter : 1.80x 0.35

Material : Titanium 5

Head Diameter : 2.45mm

Head Length : 2mm

Total Length : 9mm

Wrench type : HEX 0.50 in 1.25mm



RP & WP:

Thread Diameter : 2.00x 0.40

Material : Titanium 4

Head Diameter : 2.45mm

Head Length : 2mm

Total Length : 10.3mm

Wrench type : HEX 0.50 in 1.25mm

All Screws use the same Wrench and are 100% compatible with Nobel Replace screws.

Code	Name	Price	+
117337	XPS INTHEX Titanium Fixation Screw only for RP 3.5 0.048in (1.25mm) Hexagon	US\$ 5.85	+Info
117338	XPS INTHEX Titanium Fixation Screw only for RP 4.3 / 5.0 0.048in (1.25mm) Hexagon	US\$ 5.85	+Info

Abutments Internal Hex Ti-Base Abutments



HATHO POLISHING KIT

Thermoplastics
like *never* before



Packing Diamond D[®]
Denture Acrylic with
Tecnoflask



Keystone Industries: Laboratory




FREE SHIPPING
FOR ALL DOMESTIC ACCOUNTS!



**US
DENTAL
DEPOT**

+1 954 874 6325 
www.usdentaldepot.com 

Ti-Base Abutments - DIGITECH




Other Models

Other Models

Code	Name	Price	+
470507	DIGITECH - Screw for abutment	US\$ 6.84	+Info

Ti-Base Int-Hex Non-Rotationals - DIGITECH



Code	Name	Price	+
470606	Ti-Base UNO Hex base 4mm GH1mm	US\$ 20.43	+Info
470607	Ti-Base UNO Hex base 4mm GH2mm	US\$ 20.43	+Info
470608	Ti-Base UNO Hex base 4mm GH3mm	US\$ 20.43	+Info
470609	Ti-Base UNO Hex base 4mm GH4mm	US\$ 20.43	+Info
470610	Ti-Base UNO Hex base 6mm GH1mm	US\$ 20.43	+Info
470611	Ti-Base UNO Hex base 6mm GH2mm	US\$ 20.43	+Info
470612	Ti-Base UNO Hex base 6mm GH3mm	US\$ 20.43	+Info
470613	Ti-Base UNO Hex base 6mm GH4mm	US\$ 20.43	+Info

Ti-Base Int-Hex Rotationals - DIGITECH

Code	Name	Price	+
470614	Ti-Base UNO round base 4mm GH1mm	US\$ 20.43	+Info

470615	Ti-Base UNO round base 4mm GH2mm	US\$ 20.43	+Info
470616	Ti-Base UNO round base 4mm GH3mm	US\$ 20.43	+Info
470617	Ti-Base UNO round base 4mm GH4mm	US\$ 20.43	+Info
470618	Ti-Base UNO round base 6mm GH1mm	US\$ 20.43	+Info
470619	Ti-Base UNO round base 6mm GH2mm	US\$ 20.43	+Info
470620	Ti-Base UNO round base 6mm GH3mm	US\$ 20.43	+Info
470621	Ti-Base UNO round base 6mm GH4mm	US\$ 20.43	+Info

Abutments Internal Hex

Transfers Internal Hex




TRANSFERS INTERNAL HEX - XPS

NEEDED INSTRUMENTS

Fixation Screw 1.2mm (0.48 in)

Screwdriver 0.048in 1.2 mm for RP and WP and

Screwdriver 0.035in 0.9 mm for NP



Code	Name	Price	+
117353	XPS INTHEX Impression Transfer Closed Tray for NP & RP 2.45 Engaging - compatibl...	US\$ 14.50	+Info
117354	XPS INTHEX Impression Transfer Open Tray for RP 4.0 -(Engaging 2.45 - Platform ...	US\$ 14.50	+Info
117355	XPS INTHEX Impression Transfer Closed Tray for WP 5.0 (2.45 mm Engaging - 5.0 P...	US\$ 16.00	+Info

Abutments Internal Hex

UCLA Internal Hex



UCLA ABUTMENTS INTERNAL HEX - XPS

NEEDED INSTRUMENTS

Same instruments than used for NR

Wrench (Square engagement)

Fixation Screw 1.2 mm (0.48in)

Screwdriver 0.048in 1.2 mm for RP and WP and

Screwdriver 0.035in 0.9 mm for NP



UCLA PLASTIC ENGAGIN - XPS

Code	Name	Price	+
117327	XPS INTHEX UCLA Plastic No-Rotational for NP & RP 3.5 Diameter (Engaging 2.45 -...	US\$ 10.62	+Info
117328	XPS INTHEX UCLA Plastic No-Rotational for RP 4.0 Diameter (Engaging 2.45 - Pla...	US\$ 10.62	+Info

117329	XPS INTHEX UCLA Plastic No-Rotational for WP 5.0 D 1mm Collar Height - Compati...	US\$ 10.62	+Info
--------	---	-------------------	-----------------------

UCLA TITANIUM BASE ENGAGING - XPS

Code	Name	Price	+
117598	XPS INTHEX UCLA Titanium Base No-Rotational for NP & RP 3.5 Diameter (Engaging...	US\$ 16.00	+Info
117596	XPS INTHEX UCLA Titanium Base No-Rotational for RP 4.0 Diameter (Engaging 2.45...	US\$ 16.00	+Info
117597	XPS INTHEX UCLA Titanium Base No-Rotational for WP 5.0 D (Engaging 2.45 - 5.0 P...	US\$ 16.00	+Info

Catalog Index / Indice

Abutments Internal Hex	1
Analogs Internal Hex	3
Cad/Cam Abutments	4
Healing Caps Internal Hex	5
Screws	5
Ti-Base Abutments	6
Transfers Internal Hex	8
UCLA Internal Hex	8