

Abutments TRILOBE NobelReplace

Abutment Nobel Replace



ABUTMENTS TRILOBE NOBELREPLACE

The XP-SIMPLANT assortment of Standard Compatible Abutments offers abutment solutions for most Nobel Biocare clinical indications.

The portfolio includes all the components needed to support temporary and final, single-unit to full-arch restorations.

Compatibility

Abutments and components have been designed to work with the versatility of each Nobel Biocare implant system and its respective prosthetic connection.

- * Implant device made of Titanium 4
- * Implant Carrier / Abutment made of Titanium 2
- * Fixation Screw made of Titanium 5

NEEDED INSTRUMENTS

- * Same instruments than used for NR
- * Wrench (Square engagement)
- * Fixation Screw 1.2 mm (0.48in)
- * Screwdriver 0.048in 1.2 mm for RP and WP and
- * Screwdriver 0.035in 0.9 mm for NP



STRAIGHT TITANIUM ABUTMENTS (3.5)



5.0 Platform 4.3 engagement

XPS prosthetic abutments include their fixation screw

Code	Name	Price	+
N-29476-1	XPS TRILOBE Straight Esthetic Abutment for WP/N 5.0 Diameter (AI4351/NI) (Eng...	US\$ 14.40	+Info
N-29476-3	XPS TRILOBE Straight Esthetic Abutment for WP/N 5.0 Diameter (Engagement 4.3)-...	US\$ 9.60	+Info

5.0 Platform 5.0 engagement (Nobel)

XPS prosthetic abutments include their fixation screw

Code	Name	Price	+
------	------	-------	---

N-29476-2	XPS NOBEL Replace Compatible Straight Esthetic Abutment for WP/N 5.0 Diameter (...	US\$ 9.60	+Info
N-29476-4	XPS NOBEL Replace Compatible Straight Esthetic Abutment for WP/N 5.0 Diameter (...	US\$ 9.60	+Info

TRILOBE TRILOBE

Code	Name	Price	+
N-29376	Straight Abutment 4.3D 3mm	US\$ 9.60	+Info
N-29476	Straight Abtmt 4.3D/3.5 1mm	US\$ 12.80	+Info
N-32372	Esthetic Abutment 3.5D 1mm	US\$ 12.80	+Info
N-33399	Straight Abutment 3.5D 0mm	US\$ 9.60	+Info

ESTHETIC ANGLED ABUTMENTS



15

XPS prosthetic abutments include their fixation screw



Code	Name	Price	+
N-32380	XPS TRILOBE Angled Esthetic Abutment 15 for NP/N 3.5 & RP/N 4.3 Diameter (Enga...	US\$ 12.80	+Info
N-32381	XPS TRILOBE Angled Esthetic Abutment 15 for WP/N 5.0 Diameter (Engagement 4.3 ...	US\$ 9.60	+Info

25

XPS prosthetic abutments include their fixation screw



Code	Name	Price	+
N-32381-25	XPS TRILOBE Angled Abutment 17 for RP 4.3MM - 2mm Collar - (Compatible with Nob...	US\$ 20.25	+Info
N-32380-25	XPS TRILOBE Angled Abutment 17 for NP 3.5MM - 2mm Collar - (Compatible with Nob...	US\$ 13.50	+Info

Abutments TRILOBE NobelReplace

Analogs Trilobe NobelReplace

ANALOGS TRILOBE NOBELREPLACE

Our analogs are made of Titanium

Real shape might differ from the diagram. To check real shapes please click on each individual item.

- * Implant Carrier / Abutment made of Titanium 2
- * Fixation Screw made of Titanium 5

NEEDED INSTRUMENTS

Fixation Screw 1.2mm (0.48in)
 Screwdriver 0.048in 1.2mm for RP and WP and
 Screwdriver 0.035in 0.9mm for NP



Code	Name	Price	+
N-29498	XPS TRILOBE Titanium Implant Analog for NP/N 3.5 & RP/N 4.3 Diameter - (Engagem...	US\$ 10.80	+Info
N-29500	XPS TRILOBE Titanium Implant Analog for WP/N 5.0 Diameter (Engagement 4.3)- Co...	US\$ 10.80	+Info
N-29502	XPS NR Comp. Trilobecompatible Titanium Implant Analog WP/N 5.0 Diameter (Engag...	US\$ 12.80	+Info

Abutments TRILOBE NobelReplace

Fixation Screws NobelReplace



FIXATION SCREWS TRILOBE

We have 2 versions:

- * NP (Narrow Platform)
- * RP (Regular Platform) and WP (Wide Platform)

NP:

Thread Diameter: 1.80x 0.35
 Material: Titanium 5
 Head Diameter: 2.45mm
 Head Length: 2mm
 Total Length: 9mm
 Wrench type: HEX 0.50 in 1.25mm

RP & WP:

Thread Diameter: 2.00x 0.40
 Material: Titanium 4
 Head Diameter: 2.45mm
 Head Length: 2mm
 Total Length: 10.3mm
 Wrench type: HEX 0.50 in 1.25mm

All Screws use the same Wrench and are 100% compatible with Nobel Replace screwss.



Code	Name	Price	+
------	------	-------	---

N-29475	XPS TRILOBE Titanium Fixation Screw only for RP/N 4.3 / 5.0 0.048in (1.25mm) ...	US\$ 5.20	+Info
N-29474	XPS TRILOBE Titanium Fixation Screw only for NP/N 3.5 0.048in (1.25mm) Hexagon	US\$ 5.20	+Info

Abutments TRILOBE NobelReplace

Healing Caps NobelReplace

COVER SCREWS - NOBELREPLACE

PLEASE NOTE:

XPN prosthetic and healing components utilize a 0.99 mm (.036in) Hex tool (order below).



Code	Name	Price	+
N-29433	XPS TRILOBE Cover Screw NP/N 3.5mm - 3.3mm Platform / Internal tri-channel conne...	US\$ 5.20	+Info
N-29434	XPS TRILOBE Cover Screw RP/N 4.3D 3.3mm Platform / Internal tri-channel connecti...	US\$ 5.20	+Info
N-29435	XPS TRILOBE Cover Screw WP/N 5.0D 4.3mm Platform / Internal tri-channel connecti...	US\$ 5.20	+Info

HEALING ABUTMENTS - TRILOBE

Implant device made of Titanium 4

Implant Carrier / Abutment made of Titanium 2

Fixation Screw made of Titanium 5

NEEDED INSTRUMENTS

Same instruments than used for NR

Fixation Screw 1.2mm (0.48in)

Screwdriver 0.048in 1.2 mm for RP and WP and

Screwdriver 0.035in 0.9 mm for NP



3.5mm Platform

Code	Name	Price	+
N-29436	XPS TRILOBE Healing Abutment for NP/N 3.5 & RP/N 4.3 Diameter (Engagement 3.5)...	US\$ 7.90	+Info
N-29437	XPS TRILOBE Healing Abutment for NP/N 3.5 & RP/N 4.3 Diameter (CI3555/N) (Eng...	US\$ 7.90	+Info

4.3 Platform

Code	Name	Price	+
N-33453	XPS TRILOBE Healing Abutment for WP/N 5.0 Diameter (Engagement 4.3)- 3mm Colla...	US\$ 7.90	+Info
N-33452	XPS TRILOBE Healing Abutment for WP/N 5.0 Diameter (Engagement 4.3)- 5mm Colla...	US\$ 7.90	+Info

Abutments TRILOBE NobelReplace

Transfers NobelReplace

TRANSFERS IMPLANT INTERNAL NOBELREPLACE



Closed Tray



Code	Name	Price	+
N-33537	XPS TRILOBE Impression Transfer Closed Tray for NP/N & RP/N 3.5 & 4.3 Diameter (...)	US\$ 12.80	+Info
N-33540	XPS TRILOBE Impression Transfer Closed Tray for WP/N 5.0 Diameter (Engagement 4...	US\$ 12.80	+Info

Open Tray



Code	Name	Price	+
N-29485	XPS TRILOBE Impression Transfer Open Tray for NP/N & RP/N 3.5 & 4.3 Diameter (E...	US\$ 12.80	+Info
N-33539	XPS TRILOBE Impression Transfer Open Tray for WP/N 5.0 Diameter (Engagement 4.3...	US\$ 12.80	+Info

Abutments TRILOBE NobelReplace

Uclas NobelReplace



UCLA ABUTMENTS TRILOBE - XPS

UCLA ABUTMENTS TRILOBE - XPS

XPS prosthetic and healing components utilize a 1.25mm (.050in) Hex tool (order below).

Included Items in all our prosthetic features:

- * Fixation Screw 1.25mm (0.50in)
- * Also Additional Quantities Can Be Ordered (Sold in packs of two)
- * Nobel Biocare and Nobel Replace are Trade marks of Nobel Biocare company



UCLAS ENGAGING



Engaging Castable Plastic UCLA 1mm Collar Height

Code	Name	Price	+
N-29006-2	XPS TRILOBE UCLA Plastic no-rotational for NP/N 3.5 Diameter (Engagement 3.5)-...	US\$ 11.00	+Info
N-29010-2	XPS TRILOBE UCLA Plastic no-rotational for RP/N 4.3 Diameter (Engagement 3.5)-...	US\$ 11.00	+Info
N-29014-2	XPS TRILOBE UCLA Plastic no-rotational for WP/N 5.0 Diameter (Engagement 4.3)-...	US\$ 11.00	+Info
N-29014-3	XPS TRILOBE UCLA Plastic no-rotational for WP/N 5.0 Diameter NOBEL (Engagement ...	US\$ 11.66	+Info

Engaging Cr Cobalt Base 1mm Collar Height

Code	Name	Price	+
N-29006	XPS TRILOBE UCLA CrCo/Plastic no-rotational for NP/N 3.5 Diameter (Engagement 3...	US\$ 22.60	+Info
N-29010	XPS TRILOBE UCLA CrCo/Plastic no-rotational for WP/N 5.0 Diameter (Engagement 4...	US\$ 22.60	+Info

UCLAS NO ENGAGING

Non Engaging Castable Plastic 1mm Collar Height

ff



Code	Name	Price	+
N-29007-2	XPS TRILOBE UCLA Plastic Rotational for NP/N 3.5 Diameter (Engagement 3.5)- 1m...	US\$ 11.00	+Info
N-29011-2	XPS TRILOBE UCLA Plastic Rotational for RP/N 4.3 Diameter (Engagement 3.5)- 1m...	US\$ 11.00	+Info
N-29015-2	XPS TRILOBE UCLA Plastic Rotational for WP/N 5.0 Diameter (Engagement 4.3)- 1m...	US\$ 11.00	+Info

Non Engaging Cr Cobalt Base UCLA 1mm Collar Height



Code	Name	Price	+
------	------	-------	---

	N-29007	XPS TRILOBE UCLA CrCo/Plastic Rotational for NP/N 3.5 Diameter (Engagement 3.5 ...	US\$ 24.86	+Info
	N-29011	XPS TRILOBE UCLA CrCo/Plastic Rotational for RP/N 4.3 Diameter (Engagement 3.5 ...	US\$ 24.86	+Info

Catalog Index / Indice

Abutment Nobel Replace	1
Analogs Trilobe NobelReplace	2
Fixation Screws NobelReplace	3
Healing Caps NobelReplace	4
Transfers NobelReplace	4
Uclas NobelReplace	5