

Implants

INTHEX Internal hex Implants



XGATE IMPLANT FIXATIONS



10mm

10mm

3.3mm

3.3mm

Code	Name	Price	+
463628	IntHex Implant Fixation 3.3x10mm	US\$ 74.20	+Info

3.7mm

3.7mm

Code	Name	Price	+
463633	IntHex Implant Fixation 3.75x10mm	US\$ 74.20	+Info

4.2mm

4.2mm

Code	Name	Price	+
463638	IntHex Implant Fixation 4.2x10mm	US\$ 74.20	+Info

5.0mm

5.0mm

Code	Name	Price	+
------	------	-------	---

463643	XGATE - IntHex Implant Fixation	US\$ 74.20	+Info
11.5mm			
11.5mm			
3.3mm			
3.3mm			
Code	Name	Price	+
463629	XGATE - IntHex Implant Fixation	US\$ 74.20	+Info
3.7mm			
3.7mm			
Code	Name	Price	+
463634	IntHex Implant Fixation 11.5x3.75	US\$ 74.20	+Info
4.2mm			
4.2mm			
Code	Name	Price	+
463639	IntHex Implant Fixation 4.2x11.5	US\$ 74.20	+Info
5.0mm			
5.0mm			
Code	Name	Price	+
463644	XGATE - IntHex Implant Fixation	US\$ 74.20	+Info
13mm			
13mm			
3.3mm			
3.3mm			
Code	Name	Price	+
463630	XGATE - IntHex Implant Fixation	US\$ 74.20	+Info
3.7mm			
3.7mm			
Code	Name	Price	+
463635	IntHex Implant Fixation 3.75x13mm	US\$ 74.20	+Info
4.2mm			
4.2mm			
Code	Name	Price	+
463640	IntHex Implant Fixation 4.2x13mm	US\$ 74.20	+Info
5.0mm			
5.0mm			
Code	Name	Price	+
463645	XGATE - IntHex Implant Fixation	US\$ 74.20	+Info

16mm

16mm

3.3mm

3.3mm

Code	Name	Price	+
463631	XGATE - IntHex Implant Fixation	US\$ 74.20	+Info

3.7mm

3.7mm

Code	Name	Price	+
463636	IntHex Implant Fixation 3.75x16mm	US\$ 74.20	+Info

4.2mm

4.2mm

Code	Name	Price	+
463641	XGATE - IntHex Implant Fixation	US\$ 74.20	+Info

5.0mm

5.0mm

Code	Name	Price	+
463646	XGATE - IntHex Implant Fixation	US\$ 74.20	+Info

8mm

8mm

3.3mm

3.3mm

Code	Name	Price	+
463644	IntHex Implant Fixation 3.3x8mm	US\$ 74.20	+Info

3.7mm

3.7mm

Code	Name	Price	+
463632	IntHex Implant Fixation 3.75x8mm	US\$ 74.20	+Info

4.2mm

4.2mm

Code	Name	Price	+
463637	IntHex Implant Fixation 4.2x8	US\$ 74.20	+Info

5.0mm

5.0mm

Code	Name	Price	+
463642	XGATE - IntHex Implant Fixation	US\$ 74.20	+Info

Implants

Orthodontic Micro Implants



MICRO IMPLANT - G2

Self-Drilling Thread - Temporary Anchorage Device System

Features

- MICRO IMPLANT G2 has a unique design of micro implant is self-drilling and self-tapping. Due to the design of the conical thread, drilling is eliminated in most areas of the mouth.
- In areas of high bone density, it may be necessary to utilize the 1.1mm drill provided to penetrate the cortical plate.



Code	Name	Price	+
117590	G2 - Micro Implant G2 1.6m x 8mm - ()	US\$ 12.55	+Info
117591	G2 - Micro Implant G2 1.6m x 10mm - ()	US\$ 12.55	+Info
117592	G2 - Micro Implant G2 1.6m x 12mm - ()	US\$ 13.50	+Info
117593	G2 - Micro Implant G2 2.0m x 8mm ()	US\$ 10.26	+Info
117594	G2 - Micro Implant G2 2.0mmx 10mm ()	US\$ 6.61	+Info
117595	G2 - Micro Implant G2 2.0mmx 12mm ()	US\$ 12.69	+Info
117604	G2 - Handle for Micro Implant G2 (Screwdriver - wrench) ()	US\$ 9.33	+Info

Implants

TRILOBE conection Implants



TRILOBE IMPLANTS

XPS-TRILOBE IMPLANTS

**Compatible with Nobel Replace Select
XPN Implants are made with pure Titanium 4**

RBM Implants (Reasorbable Blast Media)

RBM, or reasorbable blast media, has been utilized for over 15 years in the dental implant industry.

* This surface treatment is designed to roughen the implant surface without leaving the residual embedded blast particles or debris in the treated substrate. To achieve the desired roughening the implant is blasted with suitable hardened particles of Hydroxyapatite (beta-Tricalcium Phosphate, Alpha TcP, TTcP and Calcium Phosphate (CPP) and then subsequently dissolved from the surface with a defined passivation treatment. Surface roughness between Ra 1.2 - 1.5 is consistently achieved

* The result is a rougher implant surface then traditional acid etch treatment providing a greater surface area for osseointegration, improved retention characteristics, increased biological fixation, and maximizing implant-to-bone contact.

Hybrid design provides to the implant fixture the advantages of a textured surface and a polished 2mm collar that helps maintaining the soft tissues



MAIN FEATURES :

- * Conical design
- * Self-Tapping
- * Micro-Thread on first millimeters

- * Compressive design for immediate load techniques
- * Polished first 2 millimeters is gentle with soft tissue
- * Platform switch
- * Double packaging similar to NR
- * Packaging Color Coded
- * 3 diameters available - 3.5 4.3 & 5.0
- * Multifunction **Implant Holder/Carrier made in Titanium 2** , that works as a Carrier, Transfer and straight abutment
- * Gamma Sterilized

Implant device made of Titanium 4
 Implant Carrier made of Titanium 2
 Fixation Screw made of Titanium 5

NEEDED BURS / DRILL

All pre existing drill can be used
 We recommend XPS TRILOBE last Conical drill

NEEDED INSTRUMENTS

Same instruments than used for NR
 Wrench (Square engagement)
 Fixation Screw 1.2mm (0.48in)
 Screwdriver 0.048in 1.2mm for RP and WP and
 Screwdriver 0.035in 0.9mm for NP

NP 3.5mm

Fixation Screw 0.9mm (0.36in)

Code	Name	Price	+
402375	XPS TRILOBE Implant Fixation 3.5mm Wide x 10mm Long (Platform Switch Engagemen...	US\$ 65.10	+Info
115829	XPS TRILOBE Implant Fixation 3.5mm Wide x 11.5mm Long (Platform Switch Engagemen...	US\$ 65.10	+Info
115830	XPS TRILOBE Implant Fixation 3.5mm Wide x 13mm Long (Platform Switch Engagemen...	US\$ 68.35	+Info

RP 4.3 mm

Fixation Screw 1.2mm (0.48in)

Code	Name	Price	+
115824	XPS TRILOBE Implant Fixation 4.3mm Wide x 8mm Long (Platform Switch Engagemen...	US\$ 65.10	+Info
115825	XPS TRILOBE Implant Fixation 4.3mm Wide x 10mm Long (Platform Switch Engagemen...	US\$ 65.10	+Info
115826	XPS TRILOBE Implant Fixation 4.3mm Wide x 11.5mm Long (Platform Switch Engagemen...	US\$ 68.35	+Info
115827	XPS TRILOBE Implant Fixation 4.3mm Wide x 13mm Long (Platform Switch Engagemen...	US\$ 68.35	+Info
115828	XPS TRILOBE Implant Fixation 4.3mm Wide x 15mm Long (Platform Switch Engagemen...	US\$ 65.10	+Info

WP 5.0mm

Fixation Screw 1.2mm (0.48in)

Code	Name	Price	+
115832	XPS TRILOBE Implant Fixation 5.0mm Wide x 8mm Long (Platform Switch Engagemen...	US\$ 65.10	+Info
115831	XPS TRILOBE Implant Fixation 5.0mm Wide x 10mm Long (Platform Switch Engagemen...	US\$ 65.10	+Info
115898	XPS TRILOBE Implant Fixation 5.0mm Wide x 11.5mm Long (Platform Switch Engagemen...	US\$ 68.35	+Info

Catalog Index / Indice

INTHEX Internal hex Implants	1
Orthodontic Micro Implants	4
TRILOBE conection Implants	4